## **Folic Acid**

# IMMULITE<sup>®</sup> 2000 Folic Acid

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of folic acid in serum, heparinized plasma or ascorbic acid-treated whole blood, as an aid in clinical diagnosis and treatment of anemia.

Catalog Number: **L2KFO2** (200 tests),  
**L2KFO6** (600 tests)

Test Code: **FOL** Color: **Orange**

CDC Analyte Identifier Code: 1907  
CDC Test System Identifier Code: 10418  
CLIA Complexity Category: Moderate

### Summary and Explanation

Folic acid (folate) and vitamin B12 are nutrients essential to hematopoiesis.<sup>1</sup> Megaloblastic anemia is almost always due to lack of one of these two vitamins.<sup>1</sup> Circulating folate levels are usually normal or elevated in vitamin B12 deficiency, but red cell folate levels are frequently low in this condition.<sup>2</sup>

Folate deficiency is commonly encountered as a result of dietary deficiency (as in alcoholism) or increased demand for this vitamin (as in pregnancy).<sup>1,3</sup> Unlike vitamin B12, folate is a heat-labile vitamin susceptible to loss by prolonged cooking. Accordingly, the prevalence of folate deficiency exhibits major demographic variations, apparently reflecting differences in dietary and culinary habits.<sup>4</sup>

Circulating folate levels, being strongly influenced by recent intake, are unreliable as an index to tissue stores.<sup>2,5,6</sup> Thus, folate levels measured in serum or plasma may be normal in the face of folate deficiency. Conversely, circulating levels may be low long before tissue stores have been exhausted.<sup>2,6</sup> Accordingly, it is important to measure red cell folate levels whenever serum or plasma levels are measured.<sup>5,7,8</sup>

### Principle of the Procedure

Competitive immunoassay.

The IMMULITE 2000 performs a 2-cycle, on-board sample treatment of patient serum, plasma or ascorbic acid-treated whole blood (for measuring red cell folate). Sample, along with ligand-labeled folic acid, is first treated with dithiothreitol (DTT), in a reaction tube containing no bead, and then with sodium hydroxide/ potassium cyanide (NaOH/KCN), in a second treatment cycle. The treated sample is transferred to a second reaction tube containing a murine anti-folate binding protein antibody-coated polystyrene bead and folate binding protein (FBP). During a 30-minute incubation, folic acid released from binding proteins in the patient sample competes with ligand-labeled folic acid for binding with FBP. The bead is washed and alkaline phosphatase labeled anti-ligand is added. During the final 30-minute incubation, the alkaline phosphatase labeled anti-ligand binds to the ligand-labeled folate that was bound to the bead during the first incubation. The unbound enzyme conjugate is removed by centrifugal wash. Substrate is added and the procedure continues as described for typical immunoassays in the Operator's Manual.

**Incubation Cycles:** 4 × 30 minutes. 2 test positions per assay: 1 sample treatment cup; 1 immunoreaction cup.

**Cycles 1 and 2:** Alkaline denaturation of endogenous binding proteins.

**Cycles 3 and 4:** Protein binding and Immunoreaction.

**Folic Acid/Vitamin B12 Panel:** To achieve maximum assay throughput, patient samples for folic acid and vitamin B12 should be run alternately, beginning with a folic acid sample. See Operator's Manual on creating a panel.

### Specimen Collection

Patient must be in a fasting state.<sup>10</sup>

For whole blood and red cell folic acid determinations, use fresh heparinized or EDTA whole blood. The patient's hematocrit must be known in order to translate whole blood into red cell folic acid results, which are expressed as nanograms per milliliter (ng/mL) of packed red blood cells. Moreover, to calculate the red cell folic acid result exactly, the

patient's serum folic acid level must also be known.

EDTA plasma should not be used as a sample type.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed serum samples are inappropriate for analysis.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Folic Acid has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

#### **Serum or Heparinized Plasma**

**Volume Required:** 50 µL serum or plasma

**Storage:** 8 hours at 2–8°C, or 6–8 weeks at –20°C.<sup>12</sup> Avoid excessive exposure to direct light.

#### **Whole Blood (Red Cell Folic Acid)**

**Volume Required:** 100 µL *fresh* heparinized or EDTA whole blood (for preparation of hemolysates).

#### **Warnings and Precautions**

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large

volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chloramphenicol, at concentrations less than 0.1 g/dL has been added as a preservative. Chloramphenicol is known to cause cancer; this disclosure is required by the state of California.

**NaOH/KCN:** Solution contains cyanide. Avoid all bodily contact.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

#### **Materials Supplied**

Components are a matched set. Labels on the inside box are needed for the assay.

##### **Folic Acid Bead Pack (L2FO12)**

With barcode. 200 beads, coated with murine monoclonal anti-folic acid binding protein. Stable at 2–8°C until expiration date.

**L2KFO2:** 1 pack. **L2KFO6:** 3 packs.

##### **Folic Acid “A” Reagent Wedge (L2FOA2)**

With barcodes. Wedges (labeled "A"), each containing 15 mL dithiothreitol (DTT) in buffer, 18.1 mL ligand-labeled folic acid in buffer, and 15 mL NaOH/KCN (NH<sub>4</sub>CN) in buffer, with preservative. Prior to opening, stable at 2–8°C until expiration date. *After opening*, stable at 2–8°C for 30 days.

**L2KFO2:** 3 wedges. **L2KFO6:** 6 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

##### **Folic Acid “D” Reagent Wedge (L2FOD2)**

With barcodes. Wedges (labeled “D”) containing 11.5 mL folate binding protein in buffer; 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to anti-ligand in buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KFO2:** 1 wedge. **L2KFO6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Folic Acid Adjustors (LFOL, LFOH)

Two vials (Low and High) of lyophilized folic acid in a human protein-based matrix, with preservative. Reconstitute each vial with **3.0 mL** sterile distilled water, and mix by gentle inversion. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KFO2:** 2 sets. **L2KFO6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

### Kit Components Supplied Separately

#### Vitamin B12/Folic Acid Sample Diluent (L2FVZ)

For running red blood cell folic acid, and for on-board dilution of high serum samples. 25 mL of concentrated (ready-to-use), vitamin B12/folic acid-free human protein-based matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2FVZ:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

**Ascorbic Acid:** Freshly prepared 0.5% solution

**CON6:** Tri-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

### Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions,

adjustment, assay and quality control procedures.

Both Reagent Wedge A (L2FOA2) and Reagent Wedge D (L2FOD2) must be loaded onto the instrument prior to assay.

#### Recommended Adjustment Interval:

4 weeks.

**Quality Control Samples:** Use controls or serum pools with at least two levels (low and high) of folic acid.

### Whole Blood and Red Cell Folic Acid Determination

**Preparation of Hemolysates:** Collect the whole blood sample and determine the hematocrit. Lyse the cells in the sample by preparing a 1-in-26 dilution of the sample (fully suspended) in *freshly* prepared 0.5% ascorbic acid solution (i.e., mix 40 µL of whole blood with 1,000 µL of 0.5% ascorbic acid). Do not use sodium ascorbate.<sup>9</sup> Vortex and let stand for 180 minutes at room temperature (15–28°C) in the dark.<sup>10</sup>

**Storage of Hemolysates:** 3 hours at 2–8°C, or 2 weeks at –20°C.<sup>10,11</sup> Avoid excessive exposure to direct light.

Following the incubation period, place the treated patient sample onto a sample rack. When ordering Folic Acid for hemolysates, you must order the test as a manual dilution (X26, to account for the ascorbic acid treatment). See the IMMULITE 2000 Operator's Manual for Dilution of samples and Manual Dilutions.

### Calculations

#### Whole Blood Results

Whole blood folic acid results will be displayed with a manual dilution factor applied by entering 26 into the Manual Dilution button located in the WORKLIST menu.

(Note: The IMMULITE 2000 accounts for the manual 1-in-26 dilution made with ascorbic acid.)

#### Red Cell Results

For an approximate measure of the packed red cell folic acid, in ng/mL, multiply the whole blood folic acid concentration (R) by 100/H, where *H* is the hematocrit *in percent*.

$$\text{Red Cell Folic Acid} \approx R \times (100 / H)$$

Strictly speaking, the serum folic acid contribution should be subtracted from the whole blood folic acid concentration before multiplying by 100/*H*. Using the patient's serum folic acid level *S*, the exact equation is:

$$\text{Red Cell Folic Acid} = (R - [S \times (100 - H) / 100]) \times (100 / H)$$

The term in the square brackets is, however, in most cases so small compared to the term *R* that it may be justifiably neglected.

### Expected Values

Because folic acid levels are strongly influenced by diet and dietary supplementation, population-based reference limits can show marked demographic differences. In the United States, for example, fortification of enriched grain products, as required by the Food and Drug Administration since 1996, has led to an estimated doubling of the mean plasma folate level among subjects not using vitamin supplements,<sup>13</sup> and a decrease in the prevalence of low folate levels, i.e. levels below 3 ng/mL (7 nmol/L).<sup>13</sup> Fortification programs can be expected to have a similar impact on whole blood and red cell folate levels, which are better measures of tissue stores.

Based on its relationship to DPC's IMMULITE Folic Acid (see Method Comparison), the assay can be expected to have the same reference limits. Note: the reference range study on which the IMMULITE Folic Acid limits are based was performed in the US prior to the FDA fortification requirements.

The reference range for whole blood folic acid was determined using 122 apparently healthy subjects in a fasting state.

Serum Folic Acid	3 – 17 ng/mL (6 – 39 nmol/L)
Whole Blood Folic Acid	112 – 508 ng/mL (254 – 1,151 nmol/L)

Laboratories should establish the appropriateness of adopting these reference ranges.

### Limitations

The folate binding protein employed in the Reagent Wedge crossreacts to some extent with the "anti-folate" drug methotrexate (MTX), a compound used in the chemotherapy of cancer. (See "Specificity" section.) Since high-dose infusions of methotrexate may result in circulating levels of 45,000 ng/mL initially and levels in the order of 2,000 ng/mL 48 hours later, the kit should not be used for folic acid measurements on patients currently receiving the drug. Moreover, great caution must be exercised when interpreting folic acid levels on patients who have recently undergone therapy with methotrexate or any other drug with a similar structure.

Folic acid results near the lower limit of normal should be interpreted with caution.<sup>10</sup>

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
ng/mL × 2.266 → nmol/L

**Reportable Range:** 1 – 24 ng/mL  
(2.3 – 54 nmol/L).

**Analytical Sensitivity:** 0.8 ng/mL  
(1.8 nmol/L).

**Intraassay Precision:** Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

**Interassay Precision:** Statistics were calculated for samples assayed in 12 different runs. (See "Interassay Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three folic acid solutions (20, 70 and 180 ng/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for folic acid. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Hemolyzed samples are unsuitable for serum folate determinations.

**Lipemia:** Presence of triglycerides in concentrations up to 5,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 9 volunteers into plain, heparinized and Becton Dickinson SST<sup>®</sup> vacutainer tubes. By linear regression:

(Heparin) = 1.05 (Serum) – 0.36 ng/mL  
r = 0.996

(SST) = 0.98 (Plain Tubes) – 0.06 ng/mL  
r = 0.997

Means:

16 ng/mL (Serum)  
16 ng/mL (Heparin)  
15 ng/mL (SST)

**Method Comparison (Serum):** The assay was compared to DPC's IMMULITE Folic Acid on 156 serum samples (Concentration range: approximately 1.8 to 20 ng/mL. See graph.) By linear regression:

(IML 2000) = 0.96 (IML) – 0.54 ng/mL  
r = 0.983

Means:

8.3 ng/mL (IMMULITE 2000)  
9.2 ng/mL (IMMULITE)

## References

- 1) Herbert V. The nutritional anemias. *Hosp Pract* 1980 Mar;15(3):65-83,87-9.
- 2) Lindenbaum J. Status of laboratory testing in the diagnosis of megaloblastic anemia. *Blood* 1983;61:624-7.
- 3) Colman N. The radioisotopic investigation of anemia. *Ligand Quarterly* 1981 Fall;4(3):24-30.
- 4) Grasbeck R. Biochemistry and clinical chemistry of vitamin B12 transport and the related diseases. *Clin Biochem* 1984;17:99-107.
- 5) Cooper BA. Folic Acid: its metabolism and utilization. *Clin Biochem* 1984;17:95-8.
- 6) McNeely MDD. Folic Acid assay. In: Kaplan LA, Pesce AJ, editors. *Clinical Chemistry*. St. Louis: CV Mosby, 1984: 1402-6.
- 7) Allen RH. Clinical role and current status of serum cobalamin (vitamin B12) assays. *Ligand Quarterly* 1981 Fall;4(3):37-44,67.
- 8) Chen I-W, et al. B12. In: Kaplan LA, Pesce AJ, editors. *Clinical Chemistry*. St. Louis: CV Mosby, 1984: 1396-1400.
- 9) Netteland B, Bakke OM. Inadequate sample preparation as a source of error in determination of erythrocyte folate by competitive binding radioassay. *Clin Chem* 1977;23:1505-6.
- 10) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: WB Saunders, 1994:2056.
- 11) Kubasik NP, Graham M, Sine HE. Storage and stability of folate and vitamin B-12 in plasma and blood samples. *Clinica Chimica Acta* 1979;95:147-149.
- 12) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:246.
- 13) Jacques PF, Selhub J, Bostom AG, et al. The effect of folic acid fortification on plasma folate and total homocysteine concentrations. *New Engl J Med* 1999;340:1449-54.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

---

## Tables and Graphs


---

### Intraassay Precision (ng/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	1.6	0.11	6.9%
2	7.0	0.29	4.1%
3	10	0.24	2.4%

**Interassay Precision (ng/mL)**

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	1.6	0.14	8.8%
2	7.4	0.42	5.7%
3	11	0.57	5.2%

**Method Comparison (Serum)**

(IML 2000) = 0.96 (IML) - 0.54 ng/mL  
 $r = 0.983$

**Linearity (ng/mL)**

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	6.68	—	—
	4 in 8	3.45	3.34	103%
	2 in 8	1.78	1.67	107%
	1 in 8	1.03	<1.00	—
2	8 in 8	7.08	—	—
	4 in 8	3.68	3.54	104%
	2 in 8	1.84	1.77	104%
	1 in 8	<1.00	<1.00	—
3	8 in 8	7.96	—	—
	4 in 8	4.14	3.98	104%
	2 in 8	2.26	1.99	114%
	1 in 8	1.13	1.00	113%
4	8 in 8	11.3	—	—
	4 in 8	5.99	5.64	106%
	2 in 8	3.14	2.82	111%
	1 in 8	1.54	1.41	109%
5	8 in 8	16.1	—	—
	4 in 8	8.34	8.05	104%
	2 in 8	4.16	4.03	103%
	1 in 8	2.07	2.01	103%
6	8 in 8	16.5	—	—
	4 in 8	8.31	8.25	101%
	2 in 8	4.30	4.13	104%
	1 in 8	2.08	2.06	101%

**Specificity (ng/mL)**

Compound <sup>1</sup>	ng/mL Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
Methotrexate	100	0.9%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH