Gastrin

IMMULITE® 2000 Gastrin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of gastrin in serum, as an aid in the diagnosis and treatment of patients with disorders associated with abnormal gastrin production.

Catalog Numbers: **L2KGA2** (200 tests)

Test Code: **GAS** Color: **Light Green**

Summary and Explanation

Gastrin is a major gastrointestinal hormone. It serves to stimulate gastric acid secretion and exists in a number of molecular forms, differing from one another in the length of the polypeptide backbone and in derivatizations of individual amino acids.¹⁰ The three principal forms — G-17, G-34 and G-14 — are named for the number of amino acids which they contain.

Immunoassays for gastrin play an essential role in the identification of Zollinger-Ellison tumors (gastrinomas). These tumors are typically, but not invariably, associated with elevated gastrin levels, gastric acid hypersecretion and peptic ulcer disease.

In fasting subjects, gastrin normally circulates at levels of less than 100 pg/mL, with some day-to-day variability. (A survey of the literature indicates that the upper reference limit is to a considerable extent method-dependent, being as high as 200 or 300 pg/mL for some RIA systems.¹⁸) Fasting gastrin levels in patients with Zollinger-Ellison syndrome are typically very high, well above the reference range for healthy individuals.

Elevated values are also encountered in other conditions. Where gastric acid secretion is impaired, for example in pernicious anemia, gastrin levels are characteristically (and appropriately) increased. Hypergastrinemia and hypersecretion of gastric acid are also encountered in the absence of pancreatic or duodenal tumors. Thus,

hypergastrinemia without gastrinoma may be found in pyloric obstruction with antrum distension, after vagotomy, in the "retained antrum" syndrome, and in some patients with ordinary peptic ulcer disease.

Since roughly half of all patients with Zollinger-Ellison tumors have fasting gastrin levels less than 500 pg/mL, the range for gastrinoma overlaps significantly with the range for other forms of hypergastrinemia. A confirmation procedure is therefore often necessary. This usually involves analysis of gastrin levels following secretin injection, calcium infusion or a test meal. The literature suggests that the secretin test is the most reliable of these follow-up procedures.^{1,6,9}

The molecular heterogeneity of gastrin has important implications for the design of gastrin immunoassays. Since some gastrinomas secrete only G-17, and others secrete only G-34, it is advantageous to use antibodies that will recognize multiple forms. Use of an overly specific immunoassay carries the risk of missing a tumor.^{4,11,19}

Principle of the Procedure

IMMULITE 2000 Gastrin is a chemiluminescent, enzyme-labeled immunometric assay based on a ligand-labeled murine monoclonal capture antibody specific for gastrin and separation by anti-ligand-coated solid phase.

The patient sample along with the ligand-labeled anti-gastrin monoclonal antibody, an alkaline phosphatase-conjugated rabbit polyclonal anti-gastrin antibody and an alkaline phosphatase-conjugated murine monoclonal anti-gastrin antibody are simultaneously incubated in the presence of the immobilized anti-ligand bead in a Reaction Tube. During the 60 minute incubation, gastrin molecules in the sample form antibody sandwich complexes which, in turn, bind to anti-ligand on the solid phase. Unbound conjugate is then removed by a centrifugal wash, after which luminogenic substrate is added and the Reaction Tube is incubated for a further five minutes.

The chemiluminescent substrate, a phosphate ester of adamantyl dioxetane, undergoes hydrolysis in the presence of alkaline phosphatase to yield an unstable intermediate. The continuous production of this intermediate results in the sustained emission of light. The bound complex — and thus also the photon output, as measured by the luminometer — is proportional to the concentration of gastrin in the sample.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

For basal serum gastrin levels, the patient must be fasting overnight, preferably 12 hours or more. Collect blood by venipuncture¹⁶ into plain tubes (without anticoagulant) noting the time of collection and separate the serum from the cells in a refrigerated centrifuge as soon as possible. Aliquot and freeze without delay.¹⁷

EDTA plasma should *not* be used in the IMMULITE 2000 Gastrin procedure.

Lipemic, icteric or grossly contaminated samples may give erroneous results. The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Gastrin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum.

Storage: Stable 4 hours at 2–8°C, or 30 days aliquotted at –20°C in a freezer that does not self-defrost.²⁰

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Gastrin Bead Pack (L2GA12)

With barcode. 200 beads, coated with anti-ligand derived from streptavidin. Stable at 2–8°C until expiration date.

L2KGA2: 1 pack.

Gastrin Reagent Wedge (L2GAA2)

With barcode. 11.5 mL ligand-labeled anti-gastrin murine monoclonal antibody, alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-gastrin antibody, and alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-gastrin, in buffer. Stable at 2–8°C until expiration date.

L2KGA2: 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Gastrin Adjustors (LGAL, LGAH)

Two vials (Low and High) of lyophilized synthetic human gastrin G-17 in a buffer matrix. Reconstitute each vial by adding 2 mL distilled or deionized water. Mix by gentle swirling or inversion. Stable after reconstitution for 30 days (aliquotted)

at –20°C in a freezer that does not self-defrost.

L2KGA2: 1 set.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For on-board dilution of high samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also available

LGACM: Bi-level Gastrin Control Module

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of gastrin.

Expected Values

A reference range study was performed in which 143 fasting samples — 49 from laboratory volunteers in good health and 94 presumed normal specimens from a large reference laboratory in the United States — were assayed by the IMMULITE 2000 Gastrin procedure. This study yielded a median value of 32 pg/mL and a range, corresponding to the nonparametric central 95% of the observations, of

13 to 115 pg/mL

for fasting samples from healthy adults.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. The assay has been standardized against the Medical Research Council's Research Standard A for gastrin II, porcine [NIBSC 66/138]. (Unless otherwise specified, all results were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factors:

pg/mL × 0.47664 → pmol/L

pg/mL × 1 → mU/L [NIBSC 66/138]

Calibration Range: up to 1,000 pg/mL [NIBSC 66/138].

Analytical Sensitivity: 5 pg/mL.

High-dose Hook Effect:

None up to 226,000 pg/mL of gastrin G-17 Type I (nonsulfated).

Precision: Samples were assayed in duplicate over the course of 10 days, four runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. Because of the molecular heterogeneity of gastrin, the linearity of this assay does not demonstrate traditional characteristics of dilutional parallelism. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with two sets of three gastrin solutions (325 pg/mL, 600 pg/mL and 1,100 pg/mL) were assayed. One set contained three G-17 Type II (sulfated) solutions, and the other set contained three G-17 Type I (nonsulfated) solutions. (See "Recovery" table for representative data.)

Specificity: The antibody recognizes predominantly G-17 with lesser recognition of G-34 and mini-gastrin. (See "Specificity" table.)

Bilirubin: May interfere with the assay, causing degradation of values at levels above 50 mg/L. (See "Bilirubin" table.)

Hemolysis: Presence of hemoglobin in concentrations up to 550 mg/dL has no effect on results, within the precision of the assay. (See "Hemolysis" table.)

Lipemia: Interferes with the assay, causing degradation of values at levels above 1,000 mg/dL. (See "Lipemia" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 30 volunteers into plain glass and plastic serum, plastic sodium heparinized and EDTA tubes and plastic gel barrier tubes (SST[®]). All tubes were from Becton Dickinson. In six sets of the matched samples, equal volumes were spiked with various concentrations of gastrin, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Gastrin procedure.

EDTA is not suitable for use. EDTA plasma can average up to 50% lower values than plain serum.

(Serum glass) =
1.03 (Serum plastic) – 0.06 pg/mL
r = 0.999

(Na Heparin plastic) =
0.995 (Serum plastic) + 0.43 pg/mL
r = 0.999

(SST plastic) = 0.978 (Serum plastic) +
4.0 pg/mL
r = 0.999

Means:

128 pg/mL (Plain Serum plastic)

132 pg/ml (Plain Serum glass)

128 pg/mL (Na Heparin plastic)

129 pg/ml (SST plastic)

Method Comparison: The assay was compared to a commercially available RIA assay for gastrin (Kit A) on 100 samples from apparently healthy volunteers. (Concentration range: approximately 10 – 305 pg/mL. See graph.) By linear regression:

(IMMULITE 2000) = 1.04 (Kit A) – 24.7 pg/mL
r = 0.926

Means:

57.0 pg/mL (IMMULITE 2000)

78.3 pg/mL (Kit A)

References

- 1) Clain JE. Diagnosis and management of gastrinoma (Zollinger-Ellison syndrome). *Mayo Clinic Proceedings* 1982;57:265-8.
- 2) Fiddian-Green RG. Hypergastrinemia – what does it mean? In: Thompson NW and Vinik AI, editors. *Endocrine Surgery Update*. New York: Grune & Stratton, 1983:219-35.
- 3) Friesen SR and Tomita T. Pseudo-Zollinger-Ellison syndrome: hypergastrinemia, hyperchlorhydria without tumor. *Annals of Surgery* 1981;194:481-93.
- 4) Hesser and Kao PC. Comparison of two assays for serum gastrin (abstract). *Clinical Chemistry* 1983;29:1161-2.
- 5) Jensen RT, et al. NIH Conference. Zollinger-Ellison syndrome: current concepts and management. *Ann Intern Med* 1983;98:59-75.
- 6) Malagelada J-R, et al. Laboratory diagnosis of gastrinoma. *Mayo Clinic Proceedings* 1982;57:211-8, 219-26.
- 7) Malagelada J-R, et al. Medical and surgical options in the management of patients with gastrinoma. *Gastroenterology* 1983;84:1524-32. See also *ibid*:1621-32.
- 8) McCarthy DM. Zollinger-Ellison syndrome. *Annual Review of Medicine* 1982;33:197-215.
- 9) McGuigan JE and Wolfe MM. Secretin injection test in the diagnosis of gastrinoma. *Gastroenterology* 1980;79:1324-31.
- 10) Rehfeld JF. Four basic characteristics of the gastrin-cholecystokinin system. *American Journal of Physiology*

1981;240:G255-66. 11) Rehfeld JF, et al. Misuse of gastrin radioimmunoassay kits. *Lancet* 1983 6 Aug;ii:338. 12) Romanus ME, et al. Comparison of four provocative tests for the diagnosis of gastrinoma. *Annals of Surgery* 1983;5:608-17. 13) Stage JG and Stadil F. The clinical diagnosis of the Zollinger-Ellison syndrome. *Scandinavian Journal of Gastroenterology* 1979;53 Suppl 14:79-91. 14) Walsh JH and Lam SK. Physiology and pathology of gastrin. *Clinics in Gastroenterology* 1980 Sept;9(3):567-91. 15) Zollinger RM, Ellison EC, et al. Primary peptic ulcerations of the jejunum associated with islet cell tumors. *Annals of Surgery* 1980;192:422-30. 16) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 17) Burtis CA, Ashwood ER, Tietz textbook of clinical chemistry. Philadelphia: Saunders, 1994: 1584. 18) Lindstedt G, et al. Analytical and clinical evaluation of a radioimmunoassay for gastrin. *Clin Chem* 1985;31,1:76-82. 19) Goetze JP, Rehfeld JF, Impact of assay epitope specificity in gastrinoma diagnosis. *Clinical Chemistry*, 2003;49,2:333-334. 20) Jacobs DS, et al., Laboratory test handbook. Cleveland: Lexi-Comp, 1996: 135.

Technical Assistance

In the United States, Contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (pg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	71	3.9	5.5%	5.0	7.0%
2	111	3.5	3.2%	6.5	5.9%
3	242	15	6.2%	16	6.6%
4	469	19	4.1%	22	4.7%
5	788	30	3.8%	77	9.8%
6	894	38	4.3%	69	7.7%

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	649	—	—
	4 in 8	329	324	102%
	2 in 8	169	162	104%
	1 in 8	88.1	81.1	109%
2	8 in 8	221	—	—
	4 in 8	109	111	98%
	2 in 8	57.4	55.3	104%
	1 in 8	32.8	27.7	118%
3	8 in 8	181	—	—
	4 in 8	93.4	90.4	103%
	2 in 8	49.3	45.2	109%
	1 in 8	27.5	22.6	122%
4	8 in 8	168	—	—
	4 in 8	81.8	84.2	97%
	2 in 8	41.1	42.1	98%
	1 in 8	22.7	21.0	108%
5	8 in 8	136	—	—
	4 in 8	76.8	68.1	113%
	2 in 8	41.1	34.0	121%
	1 in 8	21.9	17.0	129%

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	36.6	—	—
	A – G-17 II	424	360	118%
	B – G-17 II	737	635	116%
	C – G-17 II	1,365	1,135	120%
	A – G-17 I	409	360	114%
	B – G-17 I	780	635	123%
2	C – G-17 I	1,443	1,135	127%
	—	22.0	—	—
	A – G-17 II	418	346	121%
	B – G-17 II	724	621	117%
	C – G-17 II	1,350	1,121	120%
	A – G-17 I	352	346	102%
3	B – G-17 I	674	621	109%
	C – G-17 I	1,379	1,121	123%

Specificity


Compound ¹	pg/mL Added ²	Apparent Conc. ³	% Cross-reactivity ⁴
Mini-Gastrin	1,000	22.9	2.3%
Gastrin G-17 Type II (sulfated)	250	313	125%
Gastrin G-17 Type II (sulfated)	500	554	111%
Gastrin G-17 Type II (sulfated)	1,000	1,103	110%
Gastrin G-34 Type I (nonsulfated)	500	170	34%
Gastrin G-34 Type I (nonsulfated)	1,000	249	25%
Gastrin G-34 Type I (nonsulfated)	2,000	433	22%
Gastrin G-34 Type II (sulfated)	50	< 10	—
Gastrin G-34 Type II (sulfated)	500	54.9	11%
Gastrin G-34 Type II (sulfated)	5,000	391	8%
Gastrin 1-13 Type I (nonsulfated)	1,000	ND	ND
Gastrin 1-13 Type I (nonsulfated)	10,000	ND	ND
Pentagastrin	10,000	ND	ND
Pentagastrin	1,000	ND	ND
Pentagastrin	500	ND	ND
Cholecystokinin (sulfated) CCK	1,000	ND	ND
Cholecystokinin (nonsulfated)	1,000	ND	ND
Cholecystokinin CCK 10-20	1,000	9.4	0.9%
Cholecystokinin CCK 1-21	1,000	ND	ND
Caerulein	10,000	ND	ND
Caerulein	1,000	ND	ND
Caerulein	500	ND	ND
Caerulein	250	6.3	2.5%

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O Sistema de Qualidade da Diagnostic Products Corporation está registrado sob ISO 13485:2003.
EC REP DPC Biermann GmbH