

High-sensitivity Robot Arms

Tim Schröder

Article from the customer magazine Medical Solutions, October 2006

www.siemens.com/medical-magazine

SIEMENS
medical

High-sensitivity Robot Arms

Robots are commonplace with car manufacturers. The buzzing servants, made from steel and aluminum, weld metal, pack door sheeting, and lift entire carriage skeletons through the factory hall. But now these tireless assistants are faced with new and unusual tasks: Starting next year, they will position patients with previously unknown sensitivity, flexibility, and highest comfort for patient and personnel.

By Tim Schröder

Robots are part of the new Heidelberg Ion Beam Therapy Center, where highest precision is the order of the day. To optimally spare adjacent tissue, the particle beam of carbon ions or protons has to impinge on the tumor with millimeter precision.

Whenever the patient has to be positioned with great accuracy, the robot will take over: for treatment planning, for position verification, and during treatment. The robot lifts the patient tabletop from below and automatically guides it into the correct position. The flexibility of the load arm, with its six degrees of freedom, is highly advantageous. This means that the arm does more than move up and

down, forward and backward, or from left to right like a common stretcher. In addition, rotational movements are possible just as with robots in the factory. As a result, it is possible to gently position the patient in front of the beam outlet so that the particle beam reaches the target tissue at the planned entrance angle with the highest accuracy.

The Components Division of Siemens Medical Solutions in Kemnath, Germany, contracted KUKA Roboter and is developing the new highly accurate positioning and verification machine. KUKA Roboter GmbH (Ltd.) was the company of choice because it successfully

established itself years ago in the automobile industry. Headquartered in Augsburg, Germany, the company produces five basic robot models. One of these was tailor-made according to the specifications of Siemens Medical Solutions for application in particle therapy. The center in Heidelberg is using a model that weighs 1,200 kilograms and lifts patient loads up to 200 kilograms.

Unique Precision

The engineers from Siemens and KUKA Roboter were faced with the challenge of teaching their heavy machinery the necessary precision. "To meet these requirements, we had to adjust a whole series of components," explains Ralph Berke, Head of Medical Robotics at KUKA Roboter. But he will not tell which – because today the precision of these new robots is still sui generis. The steel arm positions patients even better than within a 0.5 millimeter radius of absolute position accuracy. And that is not all: the robot can subsequently correct the position in one tenth of a millimeter increments.

Berke adds that "Industrial robots are usually trimmed for repeatability. It is their task to hit the same points, day in and day out, on thousands upon thousands of car bodies. But none of them will reach the accuracy of the Siemens model implemented in Heidelberg." The performance is amazing. A machine weighing more than a ton moves toward a point with an accuracy of one tenth millimeter. Anyone who has ever tried to maneuver a car into a tight parking spot can imagine what this means. Siemens and KUKA Roboter are currently setting the gold standard with these robots, Berke emphasizes.

Usually, industrial robots and people are separated by safety fences – primarily to ensure that careless workers are not accidentally captured by the machines. Conditions in the new therapy center are, however, completely different. Here, human and machine meet eye to eye. This requires carefully planned safety systems. The positioning robot is equipped with new safety software – the Safe Robot System. The robot uses this system to monitor given room limits by itself



FROM THE LAB directly into the Heidelberg Ion Beam Therapy center: Thanks to the robot, patients can be positioned much more precisely for therapy. Top and center photos: test at Siemens; bottom photo: installation of the system at the customer site

»The system is much more flexible and precise than conventional patient-positioning tables – and comparatively cost-effective.«

Thomas Haberer, PhD,
Scientific-Technical Director,
HIT GmbH, Heidelberg, Germany

and will, for example, trigger emergency stop commands on its own. This makes the entire system faster – fractions of a second that can make all the difference.

Man and Machine

For four years now, KUKA Roboter has been demonstrating that humans and robots can work very well together – their robot-controlled roller coasters are the attraction at amusement parks. The robot carousel twirls visitors through the air at the fair. A safety system ensures the safety of the riders. The robot at the hospital will, of course, react with much greater sensitivity.

At the beginning of the treatment, the robot lowers the patient tabletop to provide the patient with easy access to the table and enable the treatment personnel to immobilize him or her comfortably. Exact positioning and immobilization of the patient are necessary to obtain maximum precision during treatment as well. Should the immobilization of the patient require more time, a shuttle system can be used. The patient tabletop is then located on a gurney. The patient is immobilized in a separate room and then moved into the treatment room. The robot arm docks on the tabletop and moves the patient to the required position, without any physical effort of the health personnel.

To align the patient exactly in front of the beam outlet unit at the previously calculated position, the treatment room at the Heidelberg facility is equipped with a second robot which is anchored in the ceiling. This second robot carries an X-ray system. Elaborate anti-collision systems in the robot control prevent the collision of the cooperating machines. Immediately prior to the treatment, the target tissue is imaged and compared with the images from the treatment plan.

Enormous Flexibility

In case of deviation, the computer calculates a shift vector. After confirmation by the oncologist, the robot aligns the patient accordingly and sets the most optimal radiation position. Due to its six degrees of freedom, the robot can maneuver patients in a highly flexible

manner. This is not possible with the tables currently used in conventional radiation therapy. The robot also demonstrates its flexibility at the site of implementation: in addition to fixed-beam treatment rooms, it can be used in treatment rooms where the beam exit is part of a rotating gantry, or to move the patient in the planning CT in exactly the same orientation as he or she will be treated afterwards.

Considering its flexibility and enormous precision, the positioning robot is part of the essential equipment and unique selling proposition of the Heidelberg system, according to Thomas Haberer, PhD, Scientific-Technical Director of HIT (Heidelberg Ion Therapy) GmbH. "The system is much more flexible and precise than conventional patient-positioning tables – and comparatively cost-effective."

In addition to seizing and moving a patient table, the robot can also handle heavy quality assurance devices and verification tools. These are objects that enable measuring and fine-tuning of the particle beam. Among them is, for example, a 200-kilogram heavy measurement basin. The robot is equipped with an adapter which can be coupled to a wide variety of accessories. It lifts the measurement basin with ease and moves it into the proper position. "Without positioning systems, the measurement basin would have to be lifted onto the table with a great deal of effort and then shifted into the right position," says Haberer.

Siemens Medical Solutions and KUKA Roboter detailed the requirements to be met by a robot working at a particle therapy system – with special attention to safety and patient friendliness. The fusion of Siemens medical knowledge and KUKA Roboter's technology leads to this highly useful innovation for the field of particle therapy; an innovation enabling precise and time-efficient treatments.

Author: Tim Schröder is a biologist and former editor of the science section of the Berliner Zeitung. He is now a freelance writer in Oldenburg, Germany, and publishes regularly in scientific journals such as Spektrum der Wissenschaft, Max Planck Forschung, and Fraunhofer Magazin.

© 2006 by Siemens AG, Berlin and Munich,
All Rights Reserved

Publisher:

Siemens AG

Medical Solutions

Henkestrasse 127, D-91052 Erlangen, Germany

Responsible for Contents:

Ulrich E. Krips

Publishing Manager:

Silke Schumann

Chief Editor:

Doris Pischitz

Editorial Team:

Timo Schickler, Katja Stöcker, Trey Stuvek

Editorial Assistance:

evolo marketing gmbh

Production:

Norbert Moser

All at: Henkestrasse 127, D-91052 Erlangen, Germany

Phone: +49-9131-84-7529

Fax: +49-9131-84-4411

e-mail: editor.medicalsolutions.med@siemens.com

Design and Editorial Consulting:

independent Medien-Design, München, Germany

in cooperation with Primafila AG, Zürich, Switzerland

Art Direction:

Horst Moser

Editorial Coordination:

Stephanie Wiesner

All at:

Widenmayerstrasse 16, D-80538 München, Germany

Printer:

MEDIA_asset_pool

Waldstrasse 18, D-91054 Erlangen, Germany

Photo Credits:

Cover, Cover Story: Johannes Krömer

Special Mobile Mammography: drehmomente.de GmbH

Krems: Toni Anzensberger

Fox Chase: Johannes Krömer

Particle Therapy: Jürgen Hinterleithner

Shared Services: Wim Klerkx

Hybrid Imaging: Imke Lass, Johannes Krömer

Prostate: Marion Stephan

Russia: Sergej Maximishin

Positioning Robot: Jürgen Hinterleithner, Günther B. Kögler

Note in accordance with § 33 Para.1 of the Federal Data Protection Law: Dispatch is made using an address file which is maintained with the aid of an automated data processing system. We remind our readers that when printed, X-ray films never disclose all the information content of the original. Artifacts in CT, MR, ultrasound and DSA images are recognizable by their typical features and are generally distinguishable from existing pathology. As referenced below, healthcare practitioners are expected to utilize their own learning, training and expertise in evaluating images.

Partial reproduction in printed form of individual contributions is permitted, provided the customary bibliographical data, such as author's name and title of the contribution as well as date and pages of *Medical Solutions*, are named. The editors request that two copies be sent to their attention. The consent of the authors and editors is required for the complete reprint of an article.

Manuscripts submitted without prior agreement as well as suggestions, proposals, and information are always welcome; they will be carefully assessed and submitted to the editorial conference for attention.

Medical Solutions on the Internet:

www.siemens.com/medical-magazine

Disclaimers

Practice of Medicine: "The information presented in this magazine is for illustration only and is not intended to be relied upon by the reader for instruction as to the practice of medicine. Any healthcare practitioner reading this information is reminded that they must use their own learning, training, and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens Medical Solutions to be used for any purpose in that regard."

Contrast Agents: "The drugs and doses mentioned herein are consistent with the approved labeling for uses and/or indications of the drug. The treating physician bears the sole responsibility for the diagnosis and treatment of patients, including drugs and doses prescribed in connection with such use. The Operating Instructions must always be strictly followed when operating the MR/CT System. The source for the technical data is the corresponding data sheets."