I M M U L I T ^E
2500

Unconjugated Estriol

DPC[®]

IMMULITE[®] 2500 Unconjugated Estriol

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2500 Analyzer — for the quantitative measurement of unconjugated (free) estriol in serum, as an aid in monitoring fetal maturity and well-being in the context of high-risk and poorly dated pregnancies.

Catalog Number: **L5KEF2** (200 tests),
L5KEF6 (600 tests)

Test Code: **FE3** Color: **Light Blue**

Summary and Explanation

Most of the estriol circulating or excreted during the third trimester of pregnancy is the joint product of fetus and placenta, originating from a precursor synthesized in the fetus by the adrenal glands, and transformed by the fetal liver and the placenta into estriol. On traversing the placenta, this is rapidly metabolized, primarily in the maternal liver, to conjugated forms: the estriol sulfates and glucuronides. As a result, "free" estriol, the unconjugated form, accounts for barely nine percent of the total estriol in circulation; the estriol sulfates, which are relatively long-lived, account for roughly half. Urinary estriol consists entirely of conjugated forms since only free estriol enters the maternal circulation while only the conjugated forms are excreted.

Normally, as the fetus develops, estriol production increases, resulting in a nearly three-fold rise in circulating estriol levels during the final trimester, and a corresponding increase in urinary levels. There is typically a surge at about the 36th week. According to the literature, free and total estriol concentrations reach approximately 15 and 250 ng/mL at term, while the urinary output climbs to approximately 45 mg/day. After 40 weeks, estriol levels gradually subside, declining by roughly 12 percent per week.⁹

There is considerable patient-to-patient variability: the reference range for a given gestational age may encompass estriol levels from 50 to 200 percent of the median for that age. Hence the *pattern*

generated by serial determinations is ordinarily of greater significance than the results of isolated measurements. Persistently low or rapidly falling estriol levels suggest fetal distress. However, because estriol concentrations are subject to diurnal and episodic variation, it is common practice to refer serum measurements to a baseline, defined for the patient as either the average or the highest of her three most recent estriol results. A drop of 40 percent or more relative to this baseline is likely to be significant.

In combination with other techniques for fetal surveillance, serial determinations have been used in the management of pregnancies complicated by diabetes,^{7,9} hypertension, prolonged gestation and uncertain dates. These clinical applications have been recently reviewed.^{7,9,10,14}

Principle of the Procedure

IMMULITE 2500 Unconjugated Estriol is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation: 60 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Samples collected into certain brands of gel barrier tubes have been shown to yield responses that are biased relative to those obtained with samples collected into plain serum tubes. See Alternate Sample Type section.

High patient samples should be *manually* diluted with an appropriate diluent.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2500 Unconjugated Estriol has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 20 µL serum.

Storage: 7 days at 2–8°C, or 6 months at –20°C.

Although serum samples have been reported to be stable at –20°C, unpublished observations suggest that a small increase, of approximately 5–10%, may occur during storage at this temperature. Therefore, it is preferable to store individual samples at 2–8°C and assay them within 1 week of collection.

It is important that all samples drawn for serial monitoring of a patient's unconjugated estriol level be measured using the same assay system, and that all be subjected to the same storage conditions. Note also that between-run sources of variability can be eliminated by assaying (or reassaying) all of the samples in one and the same run

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Unconjugated Estriol Bead Pack (L2EF12)

With barcode. 200 beads, coated with polyclonal rabbit anti-estriol. Stable at 2–8°C until expiration date.

L5KEF2: 1 pack. **L5KEF6:** 3 packs.

Unconjugated Estriol Reagent Wedge (L2EFA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to estriol in buffer, with preservative. Stable at 2–8°C until expiration date.

L5KEF2: 1 wedge. **L5KEF6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Unconjugated Estriol Adjustors (LEFL, LEFH)

Two vials (Low and High), 4 mL each, of unconjugated estriol in a processed human serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L5KEF2: 1 set. **L5KEF6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2500 Operator's Manual.

See the IMMULITE 2500 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of unconjugated estriol.

Expected Values

The reference range values were generated on the IMMULITE analyzer.

Based on the relationship of IMMULITE 2500 Unconjugated Estriol to DPC's IMMULITE 2000 Unconjugated Estriol (see Method Comparison 1), and the relationship of IMMULITE 2000 Unconjugated Estriol to IMMULITE Unconjugated Estriol (see Method Comparison 2), all three assays can be expected to have essentially the same reference ranges.

A reference range study was performed using the IMMULITE Unconjugated Estriol assay on a cross-section of 268 third trimester serum samples, with gestational ages fairly evenly distributed from 27 to 40 weeks. Tabulated below, in both mass and molar units, are the medians and central 95% ranges for each gestational week, calculated by a global parametric method. The statistics show, as expected, a clear trend towards higher values as a function of gestational age.

Unconjugated Estriol, ng/mL			
Gestational Week	n	Median	Central 95% Range
27	21	6.5	2.9 – 12.7
28	21	7.3	3.3 – 14.3
29	19	8.2	3.7 – 16.0
30	19	9.2	4.1 – 17.9
31	22	10.3	4.6 – 19.9
32	18	11.4	5.1 – 22.1
33	18	12.7	5.7 – 24.4
34	20	14.0	6.3 – 27.0
35	19	15.5	7.0 – 29.7
36	20	17.0	7.7 – >30
37	19	18.7	8.5 – >30
38	20	20.4	9.3 – >30
39	18	22.3	10.2 – >30
40	14	24.3	11.1 – >30

Unconjugated Estriol, nmol/L			
Gestational Week	n	Median	Central 95% Range
27	21	23	10 – 44
28	21	25	11 – 49
29	19	29	13 – 55
30	19	32	14 – 62
31	22	36	16 – 69
32	18	40	18 – 77
33	18	44	20 – 85
34	20	49	22 – 93
35	19	54	24 – 103
36	20	59	27 – >104
37	19	65	29 – >104
38	20	71	32 – >104
39	18	77	35 – >104
40	14	84	39 – >104

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic

antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

The IMMULITE 2500 assay employs the same reagent and bead formulation as the IMMULITE 2000. All performance characteristics are equivalent between the two platforms.

All performance data except Method Comparison 1 were generated on the IMMULITE 2000 analyzer.

Conversion Factor:

ng/mL \times 3.467 \rightarrow nmol/L

Calibration Range: 0.25 – 30 ng/mL (0.87 – 104 nmol/L).

Analytical Sensitivity: 0.1 ng/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three unconjugated estriol solutions (27, 72 and 150 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for unconjugated estriol. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 18 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of estriol, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Unconjugated Estriol procedure. By linear regression:

(Heparin) = 0.95 (Serum) – 0.03 ng/mL
r = 0.988

(EDTA) = 1.01 (Serum) – 0.15 ng/mL
r = 0.961

(SST) = 1.12 (Plain Tubes) – 0.05 ng/mL
r = 0.957

Means:

4.6 ng/mL (Serum)

4.3 ng/mL (Heparin)

4.8 ng/mL (EDTA)

5.1 ng/mL (SST Tubes)

Method Comparison 1: IMMULITE 2500 Unconjugated Estriol was compared to DPC's IMMULITE 2000 Unconjugated Estriol on 64 patient samples. (Concentration range: approximately 0.15 to 30 ng/ml. See graph 1.) By linear regression:

(IML 2500) = 0.998 (IML 2000) + 0.17 ng/mL
r = 0.996

Means:

8.42 ng/mL (IMMULITE 2500)

8.28 ng/mL (IMMULITE 2000)

Method Comparison 2: IMMULITE 2000 Unconjugated Estriol was compared to DPC's IMMULITE Unconjugated Estriol on 71 samples. (Concentration range: approximately 1.0 to 26 ng/mL. See graph 2.) By linear regression:

(IML 2000) = 1.29 (IML) + 0.38 ng/mL
r = 0.996

Means:
8.0 ng/mL (IMMULITE 2000)
5.9 ng/mL (IMMULITE)

References

1) Alexander S, et al. Renal clearance of estriol and its conjugates in normal and abnormal pregnancies. *J Clin Endocrinol Metabol* 1979;49:588-93. 2) Bagger PV, Jacobsson K, Gullberg B. Early prediction of individual total estriol curves in pregnancy. *Acta Obstet Gynecol Scand* 1981;60:187-90. 3) Bashore RA, Westlake JR. Plasma unconjugated estriol values in high risk pregnancy. *Amer J Obstet Gynecol* 1977;128:371-80. 4) Buster JE. Clinical applications of steroid assay tests of fetal-placental function. In: Abraham GE, editor. *Radioassay systems in clinical endocrinology*. New York: Marcel Dekker, 1981; 349-72. 5) Buster JE, et al. Subhourly variability of circulating third trimester maternal steroid concentrations as a source of sampling error. *J Clin Endocrinol Metabol* 1978;46:907-10. 6) Buster JE, et al. Time trend analysis of unconjugated estriol concentrations in third-trimester pregnancy. *Obstet Gynecol* 1980;56:743-47. 7) Buster JE, Freeman AG, Hobel CJ. An algorithm for determining gestational age from unconjugated estriol levels. *Obstet Gynecol* 1980;56:649-55. 8) Chard T. Normality and abnormality. In: Klopper A, editor. *Plasma hormone levels in evaluation of fetal wellbeing*. Edinburgh: Churchill Livingstone, 1976; 1-19. 9) Distler W, et al. Estriol in pregnancy. V: Unconjugated and total plasma estriol in the management of pregnant diabetic patients. *Amer J Obstet Gynecol* 1978;130:424-31. 10) Gauthier RJ, Griego BD, Goebelsmann U. Estriol in pregnancy. VII: Unconjugated plasma estriol in prolonged gestation. *Amer J Obstet Gynecol* 1981;139:382-89. 11) Gibbons JM, Nagle M. Assessing fetal health when the mother is diabetic. *Contemp Obstet Gynecol* 1980;15:115-24. 12) Goebelsmann U. The uses of estriol as a monitoring tool. *Clin Obstet Gynecol* 1979 Aug;6(2):223-44. 13) Katagiri H, et al. Estriol in pregnancy. IV: Normal concentrations, diurnal and/or episodic variations, and day-to-day changes of unconjugated and total estriol in late pregnancy plasma. *Amer J Obstet Gynecol* 1976;124:272-80. 14) Klopper A. Criteria for the selection of steroid assays in the assessment of fetoplacental function. In: Klopper A, editor. *Plasma hormone levels in evaluation of fetal wellbeing*. Edinburgh: Churchill Livingstone, 1976; 20-35. 15) Klopper A. Choice of hormone assay in the assessment of fetoplacental function. In: Fuchs F, Klopper A, editors. *Endocrinology of pregnancy*. 2nd Edition. Hagerstown, Md: Harper and Row, 1977; 350-64. 16) Klopper A. Criteria for the selection of assays in the assessment of fetoplacental function. In: Abraham GE, editor. *Radioassay systems in clinical endocrinology*. New York:

Marcel Dekker, 1981; 387-401. 17) Peterson L, Kundu N. Endocrine assessment of high risk pregnancies. *Obstet Gynecol Ann* 1980;9:169-94. 18) Scommegna A, Bieniarz J. Integrated use of physical and endocrine parameters in perinatal medicine. In: Klopper A, editor. *The human placenta*. New York: Academic Press, 1980; 181-89. 19) Tulchinsky D. Use of biochemical indices in the management of high risk obstetric patients. *Clin Perinatol* 1980 Sep; 7(2):413-21. 20) Tulchinsky D. The role of maternal and fetal hormones in pregnancy. *Contemp Obstet Gynecol* 1981 Jun;17(6):113-22. 21) Vorherr H. Placental insufficiency in relation to post-term pregnancy and fetal postmaturity. *Amer J Obstet Gynecol* 1975;123:67-103. 22) Whipple MJ, et al. Estriol in pregnancy. VI: Experience with unconjugated plasma estriol assays and antepartum fetal heart rate testing in diabetic pregnancies. *Amer J Obstet Gynecol* 1979;135:764-72. 23) DiCandio JM. Urinary estriols and total serum estriol: a correlation study. *Lab Med* 1978;9:47-8. 24) Goebelsmann U, et al. Estriol in pregnancy. *Amer J Obstet Gynecol* 1973;115:795-802. 25) Townsley JD, et al. Circadian rhythms of serum and urinary estrogens in pregnancy. *J Clin Endocrinol Metabol* 1973;36:289-95.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 310.645.8200
Fax: 310.645.9999. IMMULITE:
Tel: 973.927.2828 Fax: 973.927.4101.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	1.3	0.09	6.9%	0.16	12.3%
2	1.9	0.15	7.9%	0.23	12.1%
3	5.9	0.28	4.7%	0.43	7.3%
4	8.8	0.40	4.5%	0.50	5.7%
5	17.3	0.64	3.7%	0.98	5.7%
6	23.4	0.98	4.2%	1.38	5.9%

Linearity (ng/mL)



	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	4.9	—	—
	4 in 8	2.2	2.5	88%
	2 in 8	1.1	1.2	92%
	1 in 8	0.64	0.61	105%
2	8 in 8	7.7	—	—
	4 in 8	4.1	3.9	105%
	2 in 8	2.0	1.9	105%
	1 in 8	0.98	0.96	102%
3	8 in 8	12.2	—	—
	4 in 8	5.9	6.1	97%
	2 in 8	3.1	3.1	100%
	1 in 8	1.6	1.5	107%
4	8 in 8	19.5	—	—
	4 in 8	8.7	9.8	89%
	2 in 8	4.3	4.9	88%
	1 in 8	2.0	2.4	83%
5	8 in 8	21.9	—	—
	4 in 8	9.7	11.0	88%
	2 in 8	4.6	5.5	84%
	1 in 8	2.3	2.7	85%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH