

## **IHE Integration Statement**

| <b>Vendor</b>                 | <b>Product Name</b>              | <b>Version</b>        |
|-------------------------------|----------------------------------|-----------------------|
| SIEMENS Medical Solutions     | MAGNETOM Symphony /              | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Harmony /               | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Concerto/               | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Sonata /                | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Allegra /               | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Trio /                  | <i>syngo</i> MR 2004A |
|                               | MAGNETOM Avanto /                | <i>syngo</i> MR 2004V |
|                               | MAGNETOM C! /                    | <i>syngo</i> MR 2004C |
|                               | MAGNETOM Espree /                | <i>syngo</i> MR 2005E |
|                               | MAGNETOM Symphony, A Tim System/ | <i>syngo</i> MR 2006T |
| MAGNETOM Trio, A Tim System / | <i>syngo</i> MR 2006T            |                       |

***This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles and Actors specified below:***

| <b>Integration Profiles Implemented</b> | <b>Actors Implemented</b> |
|---|---------------------------|
| Scheduled Workflow                      | Acquisition Modality      |
|   | Image Display             |
|   | Evidence Creator          |
| Patient Information Reconciliation      | Acquisition Modality      |
| Access to Radiology Information         | Image Display             |
| Consistent Presentation of Images       | Print Composer            |

***IHE at Siemens Medical Solutions:***

<http://www.Siemens.com/IHE>

***More about products from Siemens Medical Solutions:***

<http://www.SiemensMedical.com/>

***DICOM at Siemens Medical Solutions:***

<http://www.Siemens.com/DICOM>

***Date of Statement:***

March 2006