



ACUSON X300™ Ultrasound System, Premium Edition – *Women's Imaging*

Transducers

Release 6.0

www.siemens.com/healthcare

SIEMENS



CH5-2 Transducer

Part Number:	08648086
Frequency Bandwidth:	1.4 – 5.0 MHz
Selectable 2D Frequencies:	2.5, 3.6, 5.0 MHz
Selectable THI Frequencies:	1.5, 1.8, 2.2 MHz
Selectable Color Doppler Frequencies:	2.0, 2.7 MHz
Selectable PW Doppler Frequencies:	2.2, 2.5 MHz
Number of Elements:	128
Footprint:	61.1 mm
Maximum Field of View:	66°
Maximum Display Depth:	300 mm
Biopsy Guide:	Needle Guide Kit, PN 7472314
Biopsy Guide Angles:	45° (A), 30.5° (B)
Compatible with:	ACUSON X150™ ultrasound system SONOLINE G60 S™ ultrasound system SONOLINE G40™ ultrasound system
Exam Types:	Abdomen, Early OB, Fetal Echo, OB/GYN, Peripheral Vascular, Renal

Design Attributes:

- Wide bandwidth curved array transducer
- Hanafy lens transducer technology
- User-selectable MultiHertz™ multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable



C6-2 Transducer

Part Number:	07472686
Frequency Bandwidth:	1.8 – 6.0 MHz
Selectable 2D Frequencies:	3.1, 4.0, 5.0 MHz
Selectable THI Frequencies:	2.0, 2.5 MHz
Selectable Color Doppler Frequencies:	2.5, 3.3 MHz
Selectable PW Doppler Frequencies:	2.5, 3.3 MHz
Number of Elements:	192
Footprint:	73 mm
Maximum Field of View:	90°
Maximum Display Depth:	300 mm
Biopsy Guide:	Universal Needle, Stainless, PN 5360489, Guide, PN 5260489
Biopsy Guide Angles	30°
Compatible with:	ACUSON X500™ ultrasound system, SONOLINE G60 S ultrasound system, ACUSON CV70™ cardiovascular system,
Exam Types:	Abdomen, Early OB, Fetal Echo, OB/GYN, Peripheral Vascular, Renal

Design Attributes:

- Wide bandwidth curved array transducer
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable
- Curved Vector Format



C7F2 Transducer

Part Number:	10036317
Frequency Bandwidth:	1.9 – 6.4 MHz
Selectable 2D Frequencies:	3.3, 4.2, 5.0 MHz
Selectable THI Frequencies:	2.0, 2.5 MHz
Selectable Color Doppler Frequencies:	2.5, 3.1 MHz
Selectable PW Doppler Frequencies:	2.5, 3.1 MHz
Number of Elements:	192
Footprint:	61 mm
Maximum Field of View:	91°
Maximum Display Depth:	300 mm
Biopsy Guide:	Not available
Compatible with:	ACUSON X500 ultrasound system
Exam Types:	Abdomen, Early OB, Fetal Echo, OB/GYN, Renal

Design Attributes:

- *fourSight™* 4D imaging technology
- Wide bandwidth curved array transducer
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable
- Curved Vector Format



P4-2 Transducer

Part Number:	08648045
Frequency Bandwidth:	1.3 – 4.4 MHz
Selectable 2D Frequencies:	2.0, 2.5, 2.9, 3.1, 3.6 MHz
Selectable THI Frequencies:	1.3, 1.5, 1.7, 1.8, 2.0 MHz
Selectable Color Doppler Frequencies:	2.0, 2.5, 2.7 MHz
Selectable PW Doppler Frequencies:	2.0, 2.5, 2.7 MHz
Steerable CW Doppler Frequencies:	2.0 MHz
Number of Elements:	64
Footprint:	14 mm x 17 mm
Array Length:	16 mm
Maximum Field of View:	90°
Maximum Display Depth:	300 mm
Biopsy Guide:	Not available
Compatible with:	ACUSON X150 ultrasound system SONOLINE G40 ultrasound system
Exam Types:	Abdomen, Adult Echo, OB/GYN, Renal

Design Attributes:

- Wide bandwidth phased array transducer
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable



VF8-3 Transducer

Part Number:	10038590
Frequency Bandwidth:	2.5 – 10.0 MHz
Selectable 2D Frequencies:	4.0, 5.3, 7.3, 10.0 MHz
Selectable THI Frequencies:	2.7, 3.6 MHz
Selectable Color Doppler Frequencies:	3.1, 4.0 MHz
Selectable PW Doppler Frequencies:	3.1, 4.0 MHz
2D Steering Angle:	max +/- 20°
Color Doppler Steering Angle:	max +/- 20°
PW Doppler Steering Angle:	max +/- 20°
Number of Elements:	144
Footprint:	42.2 mm
Maximum Field of View:	125 mm
Maximum Display Depth:	130 mm
Biopsy Guide:	Universal Needle Guide, Stainless, PN 5260489
Biopsy Guide Angles:	47.5°
Exam Types:	Breast, Cerebrovascular, Musculoskeletal, Orthopedic, Peripheral Vascular, Testicle, Thyroid, Venous

Design Attributes:

- Wide bandwidth linear transducer
- Virtual format imaging
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable



VF10-5 Transducer

Part Number:	08648110
Frequency Bandwidth:	4.0 – 11.4 MHz
Selectable 2D Frequencies:	6.2, 8.0, 10.0 MHz
Selectable THI Frequencies:	4.0, 4.7 MHz
Selectable Color Doppler Frequencies:	5.3, 6.2 MHz
Selectable PW Doppler Frequencies:	5.3, 6.2 MHz
2D Steering Angle:	max +/- 15°
Color Doppler Steering Angle:	max +/- 12°
PW Doppler Steering Angle:	max +/- 12°
Number of Elements:	128
Footprint:	38.4 mm
Maximum Field of View:	102 mm
Maximum Display Depth:	100 mm
Biopsy Guide:	SG-3 Needle Guide Kit, PN 7296481
Biopsy Guide Angles:	52° (A), 40° (B)
Compatible with:	ACUSON X150 ultrasound system SONOLINE G40 ultrasound system
Exam Types:	Breast, Peripheral Vascular, Thyroid

Design Attributes:

- Wide bandwidth linear transducer
- Virtual format imaging
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- 2D beam steering
- Lightweight transducer with flexible cable



VF13-5SP Transducer

Part Number:	07471506
Frequency Bandwidth:	4.3 – 13.0 MHz
Selectable 2D Frequencies:	7.3, 8.9, 11.4 MHz
Selectable THI Frequencies:	5.0, 5.7 MHz
Selectable Color Doppler Frequencies:	6.2, 7.3 MHz
Selectable PW Doppler Frequencies:	6.2, 7.3 MHz
2D Steering Angle:	max +/- 12°
Color Doppler Steering Angle:	max +/- 12°
PW Doppler Steering Angle:	max +/- 12°
Number of Elements:	128
Footprint:	13 mm x 35 mm
Array Length:	25.6 mm
Maximum Field of View:	49 mm
Maximum Display Depth:	60 mm
Biopsy Guide:	Not available
Compatible with:	ACUSON X500 ultrasound system ACUSON CV70 cardiovascular system SONOLINE G60 S ultrasound system SONOLINE G50™ ultrasound system SONOLINE Elegra™ ultrasound system
Exam Types:	Breast, Intraoperative, Musculoskeletal, Superficial Musculoskeletal, Superficial Vascular, Thyroid

Design Attributes:

- Wide bandwidth linear transducer
- Virtual format imaging
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable



VF13-5 Transducer

Part Number:	10030614
Frequency Bandwidth:	4.4 – 13.0 MHz
Selectable 2D Frequencies:	7.3, 8.9, 11.4 MHz
Selectable THI Frequencies:	4.4, 5.3 MHz
Selectable Color Doppler Frequencies:	6.2, 7.3 MHz
Selectable PW Doppler Frequencies:	6.2, 7.3 MHz
2D Steering Angle:	max +/- 12°
Color Doppler Steering Angle:	max +/- 12°
PW Doppler Steering Angle:	max +/- 12°
Number of Elements:	128
Footprint:	38.4 mm
Maximum Field of View:	61 mm
Maximum Display Depth:	60 mm
Biopsy Guide:	SG-4 Needle Guide Set, PN 7296499
Biopsy Guide Angles:	53°
Compatible with:	ACUSON X150 ultrasound system SONOLINE G40 ultrasound system
Exam Types:	Breast, Intraoperative, Musculoskeletal, Superficial Musculoskeletal, Superficial Vascular, Thyroid

Design Attributes:

- Wide bandwidth linear transducer
- Virtual format imaging
- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable



EV9-4 Transducer

Part Number:	08648037
Frequency Bandwidth:	2.9 – 9.7 MHz
Selectable 2D Frequencies:	5.0, 6.7, 8.0 MHz
Selectable THI Frequencies:	2.9, 3.3 MHz
Selectable Color Doppler Frequencies:	4.0, 5.3 MHz
Selectable PW Doppler Frequencies:	4.0, 5.3 MHz
Number of Elements:	128
Footprint:	28 mm
Maximum Field of View:	126°
Maximum Display Depth:	140 mm
Biopsy Guide:	Reusable Needle Guide set, PN 10041823
Biopsy Guide Angles:	1.7°
Compatible with:	ACUSON X150 ultrasound system SONOLINE G40 ultrasound system
Exam Types:	Early OB, OB/GYN

Design Attributes:

- Wide bandwidth endovaginal transducer
- User-selectable MultiHertz multiple frequency imaging
- Lightweight transducer with flexible cable



EV9F4 Transducer

Part Number:	10349026
Frequency Bandwidth:	2.9 – 9.0 MHz
Selectable 2D Frequencies:	4.7, 6.2, 7.3 MHz
Selectable THI Frequencies:	2.9, 3.3 MHz
Selectable Color Doppler Frequencies:	4.2, 5.3 MHz
Selectable PW Doppler Frequencies:	4.2, 5.3 MHz
Number of Elements:	192
Footprint:	38.9 mm (scan direction) x 19.9 mm (sweep direction)
Maximum Field of View:	172°
Maximum Display Depth:	140 mm
Biopsy Guide:	Needle Guide Kit, PN 10031556
Biopsy Guide Angles:	4.2°
Exam Types:	Early OB, OB/GYN

Design Attributes:

- *fourSight* 4D imaging technology
- Wide bandwidth endovaginal volume transducer
- User-selectable MultiHertz multiple frequency imaging
- Lightweight transducer with flexible cable
- Curved Vector Format



EC9-4 Transducer

Part Number:	10030334
Frequency Bandwidth:	2.9 – 9.8 MHz
Selectable 2D Frequencies:	5.0, 6.7, 8.0 MHz
Selectable THI Frequencies:	2.9, 3.3 MHz
Selectable Color Doppler Frequencies:	4.0, 5.3 MHz
Selectable PW Doppler Frequencies:	4.0, 5.3 MHz
Number of Elements:	128
Footprint:	26.2 mm
Maximum Field of View:	126°
Maximum Display Depth:	140 mm
Biopsy Guide:	Disposable Needle Guide Set, PN 10041822
Biopsy Guide Angles:	2°
Compatible with:	ACUSON X150 ultrasound system SONOLINE G40 ultrasound system
Exam Types:	Early OB, OB/GYN

Design Attributes:

- User-selectable MultiHertz multiple frequency imaging
- Ergonomically designed form factor
- Lightweight transducer with flexible cable

Frequency Bandwidth measurements represent bandwidth at +/-20 dB.

ACUSON, CV70, *fourSight*, Elegra, G40, G50, G60 S, MultiHertz, SONOLINE, X150, X300 and X500 are trademarks of Siemens Medical Solutions USA, Inc.

Local Contact Information

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355-1406 USA
Telephone: +1-888-826-9702
www.usa.siemens.com/healthcare

Europe: + 49 9131 84-0
Asia Pacific: + 65 6490 6000

Global Business Unit Address

Siemens Medical Solutions USA, Inc.
Ultrasound
1230 Shorebird Way
Mountain View, CA 94043 USA
Telephone: +1-888-826-9702
www.siemens.com/healthcare

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Munich
Germany

**Global Siemens
Healthcare Headquarters**

Global Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Telephone: + 49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
1230 Shorebird Way
Mountain View, CA 94043 USA
Telephone: +1-888-826-9702
www.siemens.com/healthcare