



Breast biopsy - efficient, improved precision,  
intuitive and patient friendly

MammoTest

Answers for life.

**SIEMENS**



## MammoTest for best access, excellent imaging and increased comfort in breast biopsy

At Siemens we know how challenging a breast biopsy can be for both you and your patient. This is why we have integrated MammoTest™ - a fully digital prone table biopsy system - into our portfolio of comprehensive breast care solutions.

Read about the benefits the Siemens MammoTest can bring to you and your facility.

### Full access helps you trap the enemy

The MammoTest breast biopsy system is designed with maximum patient and lesion access in mind. No matter how challenging or hard-to-reach the lesion seems, MammoTest helps you to literally surround the "enemy" and sample it as tissue sparing and as efficiently as possible. Throughout the procedure the patient remains in a more comfortable prone position.

### Polar coordinates with angled targeting offer best possible chest wall access

The special table design together with a rotating gantry as well as the lateral arm offer true all-angles and full 360° patient and lesion access.

### A lateral arm that comes standard

The lateral arm comes standard with every system. It provides needle access from a lateral approach, with the same accuracy and stability of the straight-on approach, which results in improved and unobstructed access to the lesion. The lateral arm makes the system especially suited for thin compressed breasts, six o'clock lesions and superficial lesions.

### Largest breast aperture in the industry

The table features a central breast aperture of 11 inches / 28 cm, which is the largest aperture available in the industry. This offers comfortable arm-through procedures ideally suited for challenging lesions close to the chest wall or in the axillary region.

### Wide range of gantry movement avoids patient repositioning

The gantry itself moves by one click only into every desired position and helps assure the shortest way for your needle. Wherever the lesion is located - you generally don't have to reposition your patient, making the entire procedure even more comfortable and efficient.

### Patient friendly from all sides

MammoTest is our digital prone table breast biopsy system. It provides best access, excellent imaging and increased comfort in breast biopsy. The system offers all the advantages of a minimally invasive intervention for a more precise and efficient diagnosis.

### Efficiency and comfort for the user

Apart from all its convincing features, MammoTest also offers compatibility with a large variety of needles and wire localisation devices. Combined with a broad range of tools and AutoGuide, the whole process of breast biopsy gets even more efficient. MammoTest is PACS and DICOM compatible, for high connectivity and increased workflow efficiency.

**For more information on the MammoTest contact [sp.med.gb@siemens.com](mailto:sp.med.gb@siemens.com)**

#### Head Office

Siemens plc  
Healthcare Sector  
Sir William Siemens Square  
Frimley, Camberley  
Surrey, GU16 8QD  
Telephone: 01276 696000  
[www.siemens.co.uk/healthcare](http://www.siemens.co.uk/healthcare)

#### Local Contact

Siemens Healthcare  
Sir William Siemens Square  
Frimley, Camberley  
Surrey, GU16 8QD  
Telephone: 01276 696000  
E-mail: [sp.med.gb@siemens.com](mailto:sp.med.gb@siemens.com)  
[www.siemens.co.uk/healthcare](http://www.siemens.co.uk/healthcare)

#### NHS Supply Chain

Emma Dunn  
Buyer - Diagnostic Imaging  
NHS Supply Chain  
Foxbridge Way  
Normanton  
West Yorkshire  
WF6 1TL  
Telephone: 01924 328728  
Mobile: 07836 264542  
Fax: 01924 328744  
E-mail: [emma.dunn@supplychain.nhs.uk](mailto:emma.dunn@supplychain.nhs.uk)  
[www.supplychain.nhs.uk](http://www.supplychain.nhs.uk)

[www.siemens.co.uk/healthcare](http://www.siemens.co.uk/healthcare)