

# IHE Integration Statement

**Product Name: ACUSON X300  
ultrasound system**

**Version: 6.0**

**Date: June 2010**

| IHE Integration Statement  |   | Date   | June, 2010 |
|--|---|--|------------|
| Vendor   | Product Name  | Version  |            |
| Siemens  | ACUSON X300™ ultrasound system  | 6.0  |            |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |   |  |            |
| Integration Profiles Implemented   | Actors Implemented  | Options Implemented  |            |
| Evidence Documents   | Evidence Creator  | OB-GYN Ultrasound Procedure Report (TID 5000)<br><br>Adult Echocardiography Procedure Report (TID 5200)<br><br>Vascular Ultrasound Procedure Report (TID 5100) |            |
| Scheduled Workflow   | Acquisition Modality  | <ul style="list-style-type: none"> <li>- Broad Worklist Query</li> <li>- Patient based Worklist Query</li> </ul>   |            |
| Echocardiography Workflow  | Acquisition Modality  | <ul style="list-style-type: none"> <li>- Broad Worklist Query</li> <li>- Patient based Worklist Query</li> </ul>   |            |
| Patient information Reconciliation   | Acquisition Modality  | <ul style="list-style-type: none"> <li>- None</li> </ul>   |            |
| <b>Internet address for vendor's IHE information:</b>  |   | <a href="http://www.siemens.com/IHE">http://www.siemens.com/IHE</a>  |            |
| <b>Links to Standards Conformance Statements for the Implementation</b>  |   |  |            |
| <b>HL7</b>   | N/A   |  |            |
| <b>DICOM</b>   | <a href="http://www.siemens.com/DICOM">http://www.siemens.com/DICOM</a> |  |            |
| <b>Links to general information on IHE</b>   |   |  |            |
| In North America:<br>www.ihe.Net   | In Europe: www.ihe-europe.org   | In Japan: www.jira-net.or.jp/ihe-j   |            |