

# IHE Integration Statement

Product Name: MAMMOMAT Inspiration

Version: AWS VB40

Date: August 2011

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare	MAMMOMAT Inspiration	AWS VB40	August 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		Assisted Acquisition Protocol Setting	
	Image Display	None	
Mammography Image	Acquisition Modality	None	
Mammography Acquisition Workflow	Acquisition Modality	No options defined	
Patient Information Reconciliation	Acquisition Modality	None	
Access to Radiology Information	Image Display	None	
Consistent Presentation of Images	Print Composer	None	
Portable Data for Imaging	Portable Media Creator	None	
	Print Composer	No options defined	
	Image Display	No options defined	
Consistent Time	Time Client	None	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	