

# IHE Integration Statement

Product Name: AXIOM Aristos  
Version: VBxx  
Date: June 2010

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	AXIOM Aristos	VBxx	June 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		Billing and Material Management
Patient Information Reconciliation	Acquisition Modality	No options defined

**Internet address for vendor's IHE information:** <http://www.siemens.com/ihe>

**Links to Standards Conformance Statements for the Implementation**

<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		

**Links to general information on IHE**

In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>
---	--	---