

IHE Integration Statement

Product Name: AXIOM Artis
Version: VB3x
Date: June 2010

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	AXIOM Artis dTA, AXIOM Artis dFA, AXIOM Artis dBA, AXIOM Artis dMP, AXIOM Artis dTC, AXIOM Artis dBC, AXIOM Artis dFC, AXIOM Artis dFC-Magnitic Navigation	VB3x	June 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		Modality Group Case
		Billing and Material Management
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query
Patient Information Reconciliation	Acquisition Modality	No options defined
Consistent Time	Time Client	None

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org
---	--	---