

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare Sector	syngo MultiModality Workplace	2008A 2008B 2008C	August 2007 April 2008 July 2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	no options defined	
	Evidence Creator	none	
Consistent Presentation of Images	Print Composer	no options defined	
Simple Image and Numeric Report	Report Creator	no options defined	
	Report Reader	no options defined	
Access to Radiology Information	Image Display	none	
Portable Data for Imaging	Portable Media Creator	none	
	Image Display	no options defined	
	Print Composer	no options defined	
	Report Reader	no options defined	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	