

# IHE Integration Statement

Product Name:  
*syngo* MultiModality Workplace

Version: 2012A  
Date: October 2011

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	syngo MultiModality Workplace	2012A	October 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	No options defined	
	Evidence Creator	None	
Consistent Presentation of Images	Print Composer	None	
Simple Image and Numeric Report	Report Creator	None	
	Report Reader	None	
Access to Radiology Information	Image Display	None	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	No options defined	
	Print Composer	No options defined	
	Report Reader	No options defined	
RT Objects	Contourer	No options defined	
	Geometric Planner	No options defined	
Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	