

IHE Integration Statement

| Vendor | Product Name | Version |
|--|---|----------------|
| SIEMENS Medical Solutions | SONOLINE Elegra | V5, V6, V7 |
| <p><i>This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles and Actors specified below:</i></p> | | |
| Integration Profiles Implemented | Actors Implemented | |
| Scheduled Workflow | Acquisition Modality | |
| Patient Information Reconciliation | Acquisition Modality | |
| Presentation of Grouped Procedures | Acquisition Modality | |
| IHE at Siemens Medical Solutions: | http://www.SiemensMedical.com/IHE | |
| More about products from Siemens Medical Solutions: | http://www.SiemensMedical.com/ | |
| DICOM at Siemens Medical Solutions: | http://www.SiemensMedical.com/DICOM | |
| Date of Statement: | November 2002 | |