

IHE Integration Statement		Date	May 2007
Vendor	Product Name	Version	
Siemens	ACUSON Antares™ ultrasound system	5.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	<ul style="list-style-type: none"> - Broad Worklist Query - Patient Based Worklist Query 	
Echocardiography Workflow	Acquisition Modality	- none	
Patient Information Reconciliation	Acquisition Modality	- none	
Evidence Documents	Evidence Creator	Echo Evidence	
Internet address for vendor's IHE information: http://www.siemens.com/IHE			
Links to Standards Conformance Statements for the Implementation			
HL7	N/A		
DICOM	http://www.siemens.com/DICOM		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	