

# IHE Integration Statement

| Vendor                | Product Name   | Version                     | Date       |
|-----------------------|--|-----------------------------|------------|
| Siemens<br>Healthcare | MAGNETOM Avanto<br>MAGNETOM Symphony, A Tim System<br>MAGNETOM Trio, A Tim System<br>MAGNETOM Espree<br>MAGNETOM ESSENZA<br>MAGNETOM Verio | syngo MR<br>VB15x,<br>VC11x | March 2008 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

| Integration Profiles Implemented   | Actors Implemented  | Options Implemented                                     |
|------------------------------------|---|---|
| Radiology Scheduled Workflow       | Acquisition Modality<br><br>Image Display<br>Evidence Creator | Broad Worklist Query<br>Billing and Material Management |
| Patient Information Reconciliation | Acquisition Modality  | Patient Based Worklist Query                            |
| Access to Radiology Information    | Image Display<br>Report Reader                                | none  |
| Simple Image and Numeric Report    | Report Creator<br>Report Reader                               | Enhanced SR<br>Enhanced SR                              |
| Consistent Presentation of Images  | Print Composer  | none  |

**Internet address for vendor's IHE information:** <http://www.siemens.com/ihe>

### Links to Standards Conformance Statements for the Implementation

|              |   |  |  |
|--------------|---|--|--|
| <b>HL7</b>   |   |  |  |
| <b>DICOM</b> | <a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a> |  |  |

### Links to general information on IHE

|   |  |   |
|---|--|---|
| In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a> | In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a> | In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a> |
|---|--|---|