

IHE Integration Statement

Product Name: *syngoMammoReport*

Version: VB41

Date: July 2011

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare	<i>syngoMammoReport</i>	VB41	July 2011

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Display	No options defined
Mammography Image	Image Display	None
	Print Composer	None
Mammography Acquisition Workflow	Image Display	No options defined
Access to Radiology Information	Image Display	None
Portable Data for Imaging	Portable Media Creator	None
	Print Composer	No options defined
	Image Display	No options defined
Consistent Time	Time Client	None

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org
---	--	---