

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Somatom Sensation Open	VB19	11/2004
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		Modality Group Case	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	No options defined in IHE	
Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>		In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>