

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Novius Radiology	27.1	9/1/2005
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Order Filler/Department System Scheduler Performed Procedure Step Manager	Billing and Material Management	
Patient Information Reconciliation	Order Filler/Department System Scheduler Performed Procedure Step Manager		
Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	