

IHE Integration Statement

**Product Name: Somatom De-
finition Flash - Definition**

Version: VA30 /VA34

Date: July 2009 / January 2010

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthcare	Somatom Definition Flash	VA30	07/2009
	Somatom Definition / Somatom Definition Flash	VA34	01/2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		Modality group case	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	No options defined in IHE	
Radiation Exposure Measurement	Acquisition Modality	none	
	Dose Info Consumer (available as a standalone tool by CARE Analytics)	Except transaction "Store dose information"	
Internet address for vendor's IHE information: http://www.siemens.com/IHE			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	