

IHE Integration Statement

Product Name: *syngo.plaza*

Version: VA20A

Date: May 2010

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare Sector	syngo.plaza	VA20A	May 2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Archive	none	
	Image Manager	none	
	Image Display	no options defined	
	PPS Manager	no options defined	
	Evidence Creator	none	
Patient Information Reconciliation (PIR)	Image Archive	no options defined	
	Image Manager	no options defined	
	PPS Manager	no options defined	
Portable Data for Imaging (PDI)	Portable Media Creator	none	
	Portable Media Importer	no options defined	
	Image Display	no options defined	
	Print Composer	no options defined	
Access to Radiology Information (ARI)	Image Archive	no options defined	
	Image Manager	no options defined	
Consistent Presentation of Images (CPI)	Image Archive	no options defined	
	Image Manager	no options defined	
	Evidence Creator	no options defined	
	Image Display	no options defined	
	Print Composer	none	
Simple Image and Numeric Report (SINR)	Report Creator	none	
Mammography Image Integration Profile	Image Archive	no options defined	
	Image Manager	no options defined	
	Image Display	none	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	none	
Consistent Time (CT)	Time Client	none	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	