

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Medical Solutions	<b>syngo® Suite</b>	V30B	Aug. 2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Department System Scheduler/Order Filler (DSS / OF) <sup>*1)</sup> Image Archive Image Manager PPS Manager	none	none
Patient Information Reconciliation (PIR)	Image Archive Image Manager PPS Manager	none	none
Reporting Workflow (RWF)	Image Archive Image Manager PPS Manager	none	none
Consistent Presentation of Images (CPI)	Image Archive Image Manager	none	none
Key Image Note (KIN)	Image Archive Image Manager	none	none
Evidence Documents (ED)	Image Archive Image Manager	none	none
Simple Image and Numeric Report (SINR)	Report Repository	none	
Access to Radiology Information (ARI)	Image Archive Image Manager Report Repository	none	Multiple Sources Multiple Sources
Basic Security <sup>*2)</sup> (SEC)	Audit Record Repository Time Server	none	none
Nuclear Medicine Imaging (NMI)	Image Archive Image Manager	none	none
Portable Data for Imaging (PDI)	Portable Media Creator <sup>*3)</sup> Image Display <sup>*3)</sup> Print Composer <sup>*3)</sup>	none	no options defined no options defined
Consistent Time (CT)	Time Server Time Client	none	none
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7 DICOM</b>	n.a. <a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	

1) For syngo® suite with syngo® Workflow MLR;

for syngo® suite with syngo® Workflow SLR: see syngo® Workflow SLR IHE Integration Statement.

2) Actors only. No full support of SEC profile.

3) Supported with a syngo® Studio Advanced Workplace.