

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	syngo® Workflow SLR	27.2	Mar. 2006

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Department System Scheduler/Order Filler (DSS / OF) PPS Manager	Billing and Material Management none
Patient Information Reconciliation (PIR)	Department System Scheduler/Order Filler (DSS / OF) PPS Manager	none none

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: <http://www.rsna.org/IHE> | In Europe: <http://www.ihe-europe.org> | In Japan: <http://www.jira-net.or.jp/ihe-j>