

Drive Productivity in Hematology


AD VIA CentraLink delivers an easy solution for true laboratory network and data management consolidation.

AD VIA CentraLink seamlessly integrates any combination of instrument analyzers with a single LIS interface. This streamlines your laboratory workflow, optimizes instrument capabilities, and improves your overall process efficiency, while minimizing the strain on your staff.

AD VIA CentraLink provides a data management and laboratory process control tool that will allow your staff to improve the consistency and quality of service that you deliver. This consolidated informatics solution will assist your laboratory staff in focusing on the delivery of timely and critical data to your customers.

Streamlining processes. Maximizing efficiency. Driving productivity.

- Review all patient results, including histograms, cytograms, instrument flags, morphology flags, QCs, and manual differentials on a single screen from a single client workstation.
- Integrated manual diff pad function allows you to perform manual differentials at any client workstation using the computer keyboard with automatic review of results.
- Customize your laboratory's protocols and utilize rules-based decision criteria to request automated slides for multiple patient types based on patient results.
- Through different levels of review, quickly make decisions on results while eliminating numerous paper printouts from multiple hematology systems.
- Manage manual slide requests from the LIS with the flexibility to suppress these requests based upon the automated patient results.
- Standardize data review with quick and easy-to-set dynamic auto-verification fields for normal ranges, delta checks, instrument flags, and QC with Westgard rules.
- Manage multispecies results from multiple hematology systems on one single screen at one single client workstation.



Link up to 16 Instrument Combinations

AD VIA CentraLink®
Networking and Data Management Solutions

Siemens Medical Solutions Diagnostics

SIEMENS
medical

ADVIA CentraLink

Capture Data in One Centralized Location

Review and Edit All on One Screen

Review patient information and test results and validate sample integrity on one easy-to-read screen. Color-coded alerts notify you in real time of results requiring your immediate attention. One view tells you where you need to validate, troubleshoot, and act on results.

Review and Edit Including Morphology Results

Diff Pad Workstation

Automated Slide Requests

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Medical Solutions Diagnostics
 511 Benedict Avenue
 Tarrytown, NY 10591-5097
 USA
 Telephone: +01-914-631-8000

www.siemens.com/diagnostics

© 2007, Siemens Medical Solutions Diagnostics
 Order No. A91DX-0700470-C2-4A00
 Printed in USA
 ADVIA, ADVIA CentraLink, and all associated marks are trademarks of Siemens Medical Solutions Diagnostics. All other trademarks and brands are the property of their respective owners.