


Real Time Access to Peer Group Data for RAPIDQC & AQC Products

QC Online/RTQC WEB Instant Access, Immediate Answers


QC Online/RTQC WEB reporting provides instant internet access to statistical analysis of your individual and peer group data on demand. Receive up-to-the-minute QC reports at any time, to monitor and improve the effectiveness of your Quality Control program.

QC Online/RTQC WEB Features and Benefits

- Web-based access allows labs to effectively manage their QC process at their convenience
- Peer data is continuously updated, eliminating QC data submission deadlines
- Significantly reduces time delays by providing real time data analysis, allowing immediate corrective action
- Easy to read graphic reports may be viewed online or printed locally for efficient QC management review



Web Statistics Report



Levey-Jennings Graph

Youden Plot

QC Online/RTQC WEB

Answers for life.

SIEMENS

Real Time QC for Siemens RAPIDComm®/RAPIDLink® System Users

RAPIDComm/RAPIDLink RTQC™ Data Converter

RAPIDComm/RAPIDLink RTQC Data Converter is a custom interface that converts QC data from the RAPIDComm/RAPIDLink Data Management System into a REAL TIME file. This allows for a quick and easy upload of data into the QC Online/RTQC WEB program.

The RAPIDComm and RAPIDLink RTQC Data Converter resides on the PC or workstation and provides an effortless process for QC data conversion. The program prepares data for electronic submission, eliminating time consuming manual transcription of the QC data.

Online training available

To enroll in the QC Online/RTQC WEB program and schedule an online training session visit: www.siemens.com/qconline.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.



RAPIDPoint® 400/405 Blood Gas Analyzers deliver central lab accuracy with Point of Care efficiency.



RAPIDLab® 1200 Blood Gas Analyzers are uniquely designed to meet your escalating throughput needs.



RAPIDComm Data Management Software gives you secure workstation access to all blood gas information – institution-wide.

Siemens Healthcare Diagnostics Inc.
1717 Deerfield Road
Deerfield, IL 60015-0778
USA

www.siemens.com/diagnostics

© 2008, Siemens Healthcare Diagnostics Inc.
Order No. A91DX-0401268-C3-4A00
Printed in USA
RAPIDComm, RAPIDLab, RAPIDLink, RAPIDPoint, RTQC and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.