

Albumin im Urin

- **Assay Format:** IMMULITE (LKHA), IMMULITE 2000 (L2KHA), IMMULITE 2500 (L5KHA)
- **Messbereich:** 2,5 - 60 µg/ml
- **Analyt. Sensitivität:** 0,5 µg/ml (IMMULITE) ; 1,0 µg/ml (IMMULITE 2000 , IMMULITE 2500)

Probanden	System	n	Median (µg/min.)	95%-Bereich (µg/min.)	95. Perc. (µg/min.)	Lit.
Gesunde	IMMULITE IMMULITE 2000 IMMULITE 2500	55	6,0	n.n. – 18,0	13,0	PILKHA-10 PIL2KHA-9 PIL5KHA-4