

C-reaktives Protein (CRP)

- **Assay Format:** IMMULITE (LKCR), IMMULITE 2000 (L2KCR), IMMULITE 2500 (L5KCR)
- **Messbereich:** bis 15 mg/dl (WHO 1st IS 85/506)
- **Analyt. Sensitivität:** 0,01 mg/dl (IMMULITE , IMMULITE 2000 ; IMMULITE 2500)
- **Umrechnungsfaktor:** mg/dl x 10 → mg/l

Probanden	System	n	Median (mg/dl)	95%-Bereich (mg/dl)	97,5. Perc. (mg/dl)	Lit.
Gesunde	IMMULITE IMMULITE 2000 IMMULITE 2500	100	0,14		1,1	PILKCR-10 PIL2KCR-10 PIL5KCR-2