

IGE , gesamt

- **Assay Format:** IMMULITE (LKIE), IMMULITE 2000 (L2KIE), IMMULITE 2500 (L5KIE)
- **Messbereich:** bis 2000 IU/ml (WHO 2nd IRP 75/502)
- **Analyt. Sensitivität:** 1,0 IU/ml (IMMULITE , IMMULITE 2000 , IMMULITE 2500)
- **Umrechnungsfaktor:** IU/ml x 2,4 → ng/ml
IU/ml x 1 → kIU/l

Probanden	System	n	Median (IU/ml)	95. Perc. (IU/ml)	97,5. Perc. (IU/ml)	Lit.
Gesunde Erwachsene	IMMULITE IMMULITE 2000 IMMULITE 2500	48	20,4	87		PILKIE-7 PIL2KIE-16 PIL5KIE-4
Kinder , 0-1 a		15	6,6	29		
Kinder , 1-2 a		34	10,1	49		
Kinder , 2-3 a		29	12,9	45		
Kinder , 3-9 a		106	14,4	52		