

LH / FSH

- **Assay Format:** IMMULITE (LKLH / LKFS) ; IMMULITE 2000 (L2KLH / L2KFS)

$$\text{LH / FSH} = \text{LH (mIU/ml)} / \text{FSH (mIU/ml)}$$

Probanden	System	n	Median	95% - Bereich	97,5. Perc.	Lit.
Frauen, Follikelphase Tag 2-5	IMMULITE	221	0,6	0,15 – 1,51	1,51	24
Frauen, Follikelphase Tag 2-9		436	0,66	0,18 – 1,64	1,64	
Frauen, Follikelphase Tag 2-11 , bis zu 5d vor LH-Peak		452	0,66	0,18 – 1,45	1,45	