

PTH-intakt

- **Assay Format:** IMMULITE (LKPP), IMMULITE 2000 (L2KPP), IMMULITE 2500 (L5KPP), IMMULITE turbo (LSKPT)
- **Messbereich:** bis 2500 pg/ml
- **Analyt. Sensitivität:** 1 pg/ml (IMMULITE) ; 0,7 pg/ml (IMMULITE 2000) ; 4 pg/ml (IMMULITE turbo)
- **Umrechnungsfaktor:** pg/ml x 0,1053 → pmol/l

Probanden	System	n	Median (pg/ml)	95% - Bereich (pg/ml)	97,5. Perc. (pg/ml)	Lit.
Gesunde Erwachsene EDTA Plasma	IMMULITE	88	38	16 - 87	87	PILKPP-13
Gesunde Erwachsene Serum		88	30	11 - 67	67	
Gesunde Erwachsene	IMMULITE turbo	174	32	8 - 74	74	PILSKPT-19
Gesunde Erwachsene		60	38	11 - 72	72	
Gesunde Erwachsene EDTA Plasma	IMMULITE 2000 IMMULITE 2500	88	38	16 - 87	87	PIL2KPP-15 PIL5KPP-5
Gesunde Erwachsene Serum	IMMULITE 2000 IMMULITE 2500	88	30	11 - 67	67	PIL2KPP-15 PIL5KPP-5