

T-Uptake

- **Assay Format:** IMMULITE (LKTU)); IMMULITE 2000 (L2KTU), IMMULITE 2500 (L5KTU)
- **Messbereich:** 10 - 70 % TU

Probanden	System	n	Median (% TU)	95% - Bereich (% TU)	97,5. Perc. (% TU)	Lit.
Gesunde Probanden	IMMULITE IMMULITE 2000 IMMULITE 2500	58	28,5	24 – 35	35	PILKTU-6 PIL2KTU-11 PIL5KTU-2
Schwangere 1. Trimester	IMMULITE	39	26,0	18 - 34	34	20
Schwangere 2.,3. Trimester		101	19,0	16 - 25	25	