

## Troponin-I, STAT

- **Assay Format:** IMMULITE 2500 STAT (L5KSTI)
- **Messbereich:** bis 100 ng/ml
- **Analyt. Sensitivität:** 0,1 ng/ml

Probanden	System	n	Median (ng/ml)	99. Perc. (ng/ml)	Lit.
Gesunde	IMMULITE 2500 STAT	255	n.n.	0,2	PL5KSTI-6
Probanden	System	n	Median (ng/ml)	Cut-off (ng/ml) 99. Percentile Gesunder	Lit.
Gesunde	IMMULITE 2500 STAT			0,2	42
Probanden	System	n	Median (ng/ml)	Cut-off (ng/ml) Bei 10% VK	Lit.
Gesunde	IMMULITE 2500 STAT			0,4	42
Probanden	System	n	Median (ng/ml)	Cut-off (ng/ml) Klinisch abgeleiteter AMI cut-off	Lit.
Gesunde	IMMULITE 2500 STAT			0,9	42