

Glucagon

- **Assay Format:** DAK RIA (KGND)
 - **Messbereich:** 25 – 500 pg/ml (WHO IS 69/194)
 - **Analyt. Sensitivität:** 13 pg/ml
- Umrechnungsfaktor:** pg/ml x 0,2872 → pmol/l

Probanden	System	n	Median (pg/ml)	95%-Bereich (pg/ml)	95. Perc. (pg/ml)	Lit.
Erwachsene	Gamma-Counter	51	112	59 - 177		PIKGND-4