

TSH

- **Assay Format:** Coat-A-Count IRMA (IKTS)
- **Messbereich:** 0,15 – 60 $\mu\text{IU/ml}$
- **Analyt. Sensitivität:** 0,03 $\mu\text{IU/ml}$

Probanden	System	n	Median ($\mu\text{IU/ml}$)	abs. Bereich ($\mu\text{IU/ml}$)	Cut-off ($\mu\text{IU/ml}$)	Lit.
Erwachsene	Gamma-Counter	443	1,3	0,3 – 5,0		PIIKTS-4
Hyperthyreose		61	n.n.	< 0,15		
Hypothyreose		77	49	5,7 - 150		