







Tim Planning Suite

1. Location on protocols in Exam Explorer
2. Set-n-Go: supports in preparing the measurement and examining large regions of the body.
3. Inline Compose: automatic composing of images from divided examination regions.
4. Position local slices on composed image
5. Coupled Graphics: Allows all graphical slice positions in all orientations within one scan to be moved at the same time.
6. Auto Coil Select: All surrounding coils are selected automatically. The maximum number of coils used is limited by the number of system channels.
7. Set-n-Go Parameters



Viewing images in Viewing Task Card



Whole Spine Imaging