

IHE Integration Statement

Vendor	Product Name	Version
SIEMENS Medical Solutions	MAGNETOM Symphony /	<i>syngo</i> MR 2004A
	MAGNETOM Harmony /	<i>syngo</i> MR 2004A
	MAGNETOM Concerto/	<i>syngo</i> MR 2004A
	MAGNETOM Sonata /	<i>syngo</i> MR 2004A
	MAGNETOM Allegra /	<i>syngo</i> MR 2004A
	MAGNETOM Trio /	<i>syngo</i> MR 2004A
	MAGNETOM Avanto /	<i>syngo</i> MR 2004V
	MAGNETOM C! /	<i>syngo</i> MR 2004C
	MAGNETOM Espree /	<i>syngo</i> MR 2005E
	MAGNETOM Symphony, A Tim System/	<i>syngo</i> MR 2006T
MAGNETOM Trio, A Tim System /	<i>syngo</i> MR 2006T	

This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles and Actors specified below:

Integration Profiles Implemented	Actors Implemented
Scheduled Workflow	Acquisition Modality
	Image Display
	Evidence Creator
Patient Information Reconciliation	Acquisition Modality
Access to Radiology Information	Image Display
Consistent Presentation of Images	Print Composer

IHE at Siemens Medical Solutions:

<http://www.Siemens.com/IHE>

More about products from Siemens Medical Solutions:

<http://www.SiemensMedical.com/>

DICOM at Siemens Medical Solutions:

<http://www.Siemens.com/DICOM>

Date of Statement:

March 2006