

## **IHE Integration Statement**

| <b>Vendor</b>             | <b>Product Name</b> | <b>Version</b>        |
|---------------------------|---------------------|-----------------------|
| SIEMENS Medical Solutions | MAGNETOM Symphony / | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Harmony /  | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Concerto/  | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Sonata /   | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Allegro /  | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Trio /     | <i>syngo</i> MR 2004A |
|                           | MAGNETOM Avanto /   | <i>syngo</i> MR 2004V |
|                           | MAGNETOM C! /       | <i>syngo</i> MR 2004C |

***This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles and Actors specified below:***

| <b>Integration Profiles Implemented</b> | <b>Actors Implemented</b> |
|---|---------------------------|
| Scheduled Workflow                      | Acquisition Modality      |
|   | Image Display             |
|   | Evidence Creator          |
| Patient Information Reconciliation      | Acquisition Modality      |
| Access to Radiology Information         | Image Display             |
| Consistent Presentation of Images       | Print Composer            |

***IHE at Siemens Medical Solutions:***

<http://www.Siemens.com/IHE>

***More about products from Siemens Medical Solutions:***

<http://www.SiemensMedical.com/>

***DICOM at Siemens Medical Solutions:***

<http://www.Siemens.com/DICOM>

***Date of Statement:***

August 2004