

Siemens		SignalAD						CAS Innovations		Medtronic		Pirion		Stryker	
C-Arm Model	Software Version	Spiral 19-24	Fluoro 14 14	Fluoro 20 20	Fluoro 24 24	Fluoro 28 28	ACL 14	Hp 14	Application						
									Fluoro 14	Fluoro 20	Fluoro 24	Fluoro 28			
<b>3D Fluoro, analog</b>															
SIEMENS Compact L1 7"	n.a.							n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
SIEMENS Compact L1 7"	n.a.							n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
SIEMENS Iso-C LF, 11, 11"	n.a.							n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
SIEMENS Compact L1 9"	n.a.	x	x	x	x	x	x	x	Fluorobase 2.2.2 and higher				Fluorobase Software 1.0 and higher (P only)		
SIEMENS Compact L1 9"	n.a.	x	x	x	x	x	x	x	Fluorobase 2.2.2 and higher	Fluorobase Advanced			Fluorobase Software 1.0 and higher (P only)		
SIEMENS Iso-C, 11, 9", 11"	n.a.	x	x	x	x	x	x	x	Fluorobase 2.2.2 and higher	Fluorobase Advanced			Fluorobase Software 1.0 and higher		
SIEMENS Iso-C/D, 11, 9" (2D Fluoro mode)	n.a.	x	x	x	x	x	x	x	Fluorobase 2.2.2 and higher	Fluorobase Advanced			Fluorobase Software 1.0 and higher		
SIEMENS Iso-C, 11, 9"	n.a.	x	x	x	x	x	x	x	Fluorobase 2.2.2 and higher	Fluorobase Advanced			Fluorobase Software 1.0 and higher		
<b>3D Fluoro, digital - NavisLink 2D</b>															
ARCADIS Vaski, 11, 9"	1010 1011 1012C	x	x	x	x	x	x	x	Fluorobase M02 1.2 and higher Spiralbase 1.4 and higher		n.a.		Fluorobase Software 2.0 and higher		
ARCADIS Orbis, 11, 9"	1010 1011 1012C 1013A	x	x	x	x	x	x	x	Fluorobase M02 1.2 and higher Spiralbase 1.4 and higher		n.a.		Fluorobase Software 2.0 and higher		
<b>3D Fluoro, digital - NavisLink 3D</b>															
SIEMENS Iso-C 3D, 11, 9" (in 3D mode)	1030B 1030D		x	x					CAPPA-System 3D CAPPA-C-View 1.0	Base 1.1 and higher Spiralbase 1.1 and higher	Fluorobase 3D		Spiral 1.1 Iso-C Software and higher		
ARCADIS Orbis 3D, 11, 9" (in 3D mode)	1010 1011 1012C 1013A								CAPPA-System 3D CAPPA-C-View 1.0	Base 1.2 and higher Spiralbase 1.2 and higher	Fluorobase 3D		Spiral 1.2 Iso-C Software and higher		