

Data



MODULARIS

Modular design with new flexibility

MODULARIS

Accessories



Accessories (optional)

Coxafix leg supports

44 98 445

The leg supports can be attached anywhere along the lateral OR rails with a quick-action lock. Continuous, gas spring assisted height adjustment with self-locking spindle, padded leg rests, lockable and swivelable in 3 directions for pressure-free positioning of the legs in the knee joint. Leg rest quick-action lock, suitable for use of the special pediatric leg rest.

Pediatric leg supports (Coxafix)

44 16 504

For examinations of children, the adult leg rests of the Coxafix leg supports can be replaced by pediatric leg rests.

Amatech foot supports

Amatech adjustable gas-spring loaded foot supports. Attached with low profile siderail clamps.

Boot stirrups, Siemens customization, including clamps and boot pads.

Adjust at side rail socket from low to high lithotomy.

The telescoping bar permits simple length adjustment while eliminating any protrusions beyond the boot. An offset support creates unrestricted access for the operative team. A limited motion ball socket at the base of the boot prevents unsafe inversion or eversion when the socket is loosened.

It may be left loosened to self align during position changes and then locked in final position.

Amatech knee crutch stirrups

For easy attachment to the table rails.

The unique "Lift Assist" neutralizes leg weight, allowing easy and quick adjustment of the leg position.

Knee crutch stirrups, Siemens customization, including clamps and crutch pads. Attachment to the lateral OR rails.

Provide excellent patient comfort and precise surgical positioning for short duration procedures. Doughty soft pad reduces popliteal pressure.

One easy-to-use handle provides vertical adjustment from 40 cm (16") to 52 cm (21") and 360° horizontal rotation.

Velcro strap for security.

Accessories (optional)

Patient arm rest (rigid)

44 98 452

The padded arm rest allows safe and comfortable positioning of the patient's left or right forearm. It can be attached and locked in place anywhere along the lateral OR rails.

Patient arm rest (flexible)

144 04 724

The padded arm rest allows safe and comfortable positioning of the patient's left or right forearm. It can be attached and locked in place anywhere along the lateral OR rails.

Shoulder supports

44 16 546

The shoulder supports allow safe and comfortable positioning of the patient in Trendelenburg positions. The padding can be adapted to the patient's anatomy. The shoulder supports attach anywhere along the lateral OR rails.

Wrist restraints

44 16 587

The wrist restraints are attached to the OR rails and are used to restrain the patient's hands during an examination.

Elbow supports

144 00 881

Swiveling elbow supports at the perineal end of the tabletop, supporting the examiner's elbows during transurethral procedures.



MODULARIS

Accessories



Accessories (optional)

Foot switch

44 98 486

Multifunctional foot control unit for all patient table movements, 4 pedals, with selector switch for 8 control functions.

Separable lower body radiation protection

44 16 512

To protect the examiner's lower body from scattered radiation, a radiation protection shield with 0.5 mm lead equivalent can be attached. Adapter slides onto accessory rail on left/right.

Transmobile

44 24 540

This upgrade kit comprises all the parts necessary for optimizing the MODULARIS Uro patient table for use in several hospitals. The larger wheels, separate brakes and a floor plate to protect the table make transport considerably easier.

Accessories (optional)

Holder for plastic drain bag

44 98 478

Flexible, tiltable spring band frame for use at the perineal table end, allowing easy attachment of a plastic drain bag with drain hose.

Urological PE drain bag (box of 20)

44 16 579

Preassembled plastic drain bags with apron, removable filter screen and expandable drain hose.



Positioning aid for infants⁶⁾

44 16 595

Comfortable positioning aid for ESWL treatment of infants. Easy-to-clean and mechanically stable plastic cradle for safe and comfortable positioning of infants during lithotripsy treatment. The positioning aid is suitable for children aged between 0 to 9 years with a height of approx. 50 cm to 110 cm.

⁶⁾ Not available in the U.S.

MODULARIS

Accessories



Accessories (optional)

Anesthesia screen holder

44 16 520

The telescopic height-adjustable anesthesia screen holder allows the attachment of a sterile cover and mounts anywhere along the OR rails.

Infusion bottle holder

144 04 722

The telescopic extendable infusion bottle holder can hold up to 4 infusion bottles and attaches anywhere along the OR rails of the table.

Paper roll holder

44 98 460

Attaches to the examination table for holding medical crepe paper.

Medical crepe paper

44 16 553

Medical crepe paper is placed between the patient and the examination table.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

The required features should therefore be specified in each individual case at completion of contract.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Original images always lose a certain amount of detail when reproduced.

This product does not yet have a certificate of conformity.

This product is not commercially available in the United States (U.S.A.).

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources, waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for new components, we guarantee the quality of these recycled components.

Siemens AG
Wittelsbacherplatz 2
D-80333 Muenchen
Germany

Headquarters

Siemens AG, Medical Solutions
Henkestr. 127, D-91052 Erlangen
Germany
Telephone ++49 9131 84-0
www.siemens.com/medical

Contact addresses:

In the USA

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355
Telephone: +01 610 448 4500
Telefax: +01 610 448 1620

In Germany

Siemens AG, Medical Solutions
Special Systems
Allee am Röthelheimpark 2
D-91052 Erlangen
Germany
Telephone ++49 9131 84-0
siemens.com/medical

Siemens **Medical**
Solutions that help