

IHE Integration Statement

**Product Name: ACUSON S2000
ultrasound system**

Version: VB10

Date: June 2010

IHE Integration Statement		Date	June 2010
Vendor	Product Name	Version	
Siemens	ACUSON S2000™ ultrasound system	VB10 (2.0)	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Evidence Documents	Evidence Creator	OB-GYN Ultrasound Procedure Report (TID 5000) Adult Echocardiography Procedure Report (TID 5200) Vascular Ultrasound Procedure Report (TID 5100)	
Scheduled Workflow	Acquisition Modality	<ul style="list-style-type: none"> - Broad Worklist Query - Patient based Worklist Query 	
Patient information Reconciliation	Acquisition Modality	<ul style="list-style-type: none"> - None 	
Internet address for vendor's IHE information: http://www.siemens.com/IHE			
Links to Standards Conformance Statements for the Implementation			
HL7	N/A		
DICOM	http://www.siemens.com/DICOM		
Links to general information on IHE			
In North America: www.ihe.Net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	