

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	ARCADIS Varic / Orbic / Orbic 3D ARCADIS Avantic	VB11 VB13B	June 2005
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
Scheduled Workflow	Image Display	none	
Patient Information Reconciliation	Acquisition Modality	none	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	none	
Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	