

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Mammomat Novation DR	AWS VA11	June 2005

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query
Scheduled Workflow	Image Display	none
Patient Information Reconciliation	Acquisition Modality	none
Access to Radiology Information	Image Display	none
Consistent Presentation of Images	Print Composer	none

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j
---	--	---