

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Mammomat Novation DR	AWS VA20	September 2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
	Image Display	None	
Mammography Image	Acquisition Modality	None	
Patient Information Reconciliation	Acquisition Modality	None	
Access to Radiology Information	Image Display	None	
Consistent Presentation of Images	Print Composer	None	
Portable Data for Imaging	Portable Media Creator	None	
	Print Composer	no options defined	
	Image Display	no options defined	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	