

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	AXIOM Artis TA, AXIOM Artis FA, AXIOM Artis BA, AXIOM Artis MP AXIOM Artis dTA, AXIOM Artis dFA, AXIOM Artis dBA, AXIOM Artis dMP AXIOM Artis TC, AXIOM Artis FC, AXIOM Artis BC AXIOM Artis dTC, AXIOM Artis dBC AXIOM Artis dFC, AXIOM Artis dFC-Magnetic Navigation	VB22	December 2005

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol setting Modality Group Case Billing and Material Management
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query
Patient Information Reconciliation	Acquisition Modality	none

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j
---------------------------------------------------------------------------------	------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------