

IHE Integration Statement

| Vendor | Product Name | Version | Date |
|--|--|--|---------------|
| Siemens Medical Solutions | AXIOM Iconos R200 (Fluorospot Compact) | VD46 | December 2005 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol setting Billing and Material Management | |
| Patient Information Reconciliation | Acquisition Modality | none | |
| Internet address for vendor's IHE information: http://www.siemens.com/ihe | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| HL7 | n.a. | | |
| DICOM | http://www.siemens.com/dicom | | |
| Links to general information on IHE | | | |
| In North America: http://www.rsna.org/IHE | In Europe: http://www.ihe-europe.org | In Japan: http://www.jira-net.or.jp/ihe-j | |