

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Artis zee / zeego	VC13	July 2008

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		Assisted Acquisition Protocol setting
		Modality Group Case
		Billing and Material Management
	Image Display	no options defined
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
	Image Display	no options defined
Patient Information Reconciliation	Acquisition Modality	no options defined
Access to Radiology	Image Display	Multiple Sources
Consistent Presentation of Images	Print Composer	no options defined
Consistent Time	Time Client	SNTP

**Internet address for vendor's IHE information:** <http://www.siemens.com/ihe>

### Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		

### Links to general information on IHE

In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>
---	--	---