

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	Artis zee / zeego	VC14	July 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		Assisted Acquisition Protocol setting	
		Modality Group Case	
		Billing and Material Management	
	Image Display	No options defined	
	Evidence Creator	None	
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
	Image Display	No options defined	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Access to Radiology	Image Display	Multiple Sources	
Consistent Presentation of Images	Print Composer	No options defined	
Consistent Time	Time Client	SNTP	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	No options defined	
	Print Composer	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n. a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	