

IHE Integration Statement

Product Name: Artis zee / zeego

Version: VC20

Date: January 2012

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	Artis zee / zeego	VC20	January 2012

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		Modality Group Case
		Billing and Material Management
	Image Display	No options defined
	Evidence Creator	None
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
	Image Display	No options defined
Patient Information Reconciliation	Acquisition Modality	No options defined
Access to Radiology Information	Image Display	Multiple Sources
Consistent Presentation of Images	Print Composer	None
Consistent Time	Time Client	SNTP
Portable Data for Imaging	Portable Media Creator	None
	Image Display	No options defined
	Print Composer	No options defined
Radiation Exposure Monitoring	Acquisition Modality	No options defined

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org
---	--	---