

# IHE Integration Statement

Product Name: syngo MR

Date: July 2010

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	MAGNETOM Symphony MAGNETOM Sonata MAGNETOM Symphony/Sonata Vision MAGNETOM Trio MAGNETOM C! MAGNETOM Harmony MAGNETOM Harmony Expert MAGNETOM Concerto MAGNETOM Allegra MAGNETOM Novus 1.5 MAGNETOM ESSENZA MAGNETOM Avanto MAGNETOM Symphony, A TIM System MAGNETOM Trio, A TIM System MAGNETOM Espree MAGNETOM Verio MAGNETOM Aera MAGNETOM Skyra	Syngo MR VA35x, VB17x, VC15x VD11x	July 2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality	Broad Worklist Query Billing and Material Management	
	Image Display Evidence Creator		
Patient Information Reconciliation	Acquisition Modality	Patient Based Worklist Query	
Access to Radiology Information	Image Display Report Reader	none	
Simple Image and Numeric Report	Report Creator Report Reader	Enhanced SR Enhanced SR	
Consistent Presentation of Images	Print Composer	none	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	