

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Magnetom Symphony, A Tim System Magnetom Trio, A Tim System Magnetom Avanto Magnetom Espree	Syngo VB13	August 2006

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality Image Display Evidence Creator	Broad Worklist Query Billing and Material Management
Patient Information Reconciliation	Acquisition Modality	Patient Based Worklist Query
Access to Radiology Information	Image Display Report Reader	none
Simple Image and Numeric Report	Report Creator Report Reader	Enhanced SR Enhanced SR
Consistent Presentation of Images	Print Composer	none

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7			
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j
---	--	---